DEVILBISS® PULSEDOSE® SAUERSTOFF-SPARSYSTEM - SERVICEHANDBUCH

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DeVilbiss Healthcare Pty. Limited
15 Carrington Road, Unit 8
Castle Hill NSW 2154
AUSTRALIA
+61-2-9899-3144

DeVilbiss Healthcare SAS
13/17, Rue Joseph Priestley
37100 Tours
FRANCE
+33 (0) 2 47 42 99 42

DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim
GERMANY
+49 (0) 621-178-98-0

DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA
800-338-1988 • 814-443-4881 • www.DeVilbissHealthcare.com

DeVilbiss® and PulseDose® are registered trademarks of DeVilbiss Healthcare.

© 2015 DeVilbiss Healthcare LLC. 05.15 All Rights Reserved.

LT-2029-DE Rev. A

DEVILBISS® PULSEDOSE® OXYGEN CONSERVING DEVICE SERVICE MANUAL

FOR QUALIFIED SERVICE PERSONNEL ONLY

MODEL PD1000A

CAUTION-Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

TABLE OF CONTENTS

GENERAL INFORMATION	2
Initial Inspection	2
Maintaining the PD1000A	2
General Description	2
THEORY OF OPERATION	2
INSTALLATION AND OPERATION	3
Installation	3
Operation	3
MAINTENANCE PROCEDURES	3
Preventive Maintenance	3
Cleaning	4
Cylinder Recertification Testing	4
TROUBLESHOOTING	4
SERVICE INSTRUCTIONS	5
Battery Door/ Battery Removal	5
Cover/Gauge Removal	5
Cannula Fitting Removal	5
Solenoid/Vacuum Switch Testing and Removal	5
PARTS AND ILLUSTRATIONS	6
PNEUMATIC DIAGRAM	6
UNIT SPECIFICATIONS	7
PARTS RETURN AND ORDERING POLICY	7
PROVIDER'S NOTES	8
GUIDANCE AND MANUFACTURER'S DECLARATION	9

GENERAL INFORMATION

Initial Inspection

An initial inspection should be performed on the PD1000A as soon as possible after receipt. When removed from carton, an inspection should be made for any damage due to shipping. If shipping damage has occurred, call DeVilbiss Healthcare at 800-338-1988 (814-443-4881) for replacement instructions.

Maintaining The PD1000A

The PD1000A should be periodically maintained according to the guidelines set forth in this manual. Maintenance and testing should only be done by qualified service personnel. Failure to follow the procedures set forth in this manual may void the warranty.

General Description

The PD1000A PulseDose system delivers a pulse or "bolus" of oxygen at the leading edge of inspiration. This bolus is delivered at both the proper flow and volume so that it is delivered deep into the lungs where gas exchange takes place. The PD1000A PulseDose is rate responsive from 6 BPM to 40 BPM.

THEORY OF OPERATION

The PD1000A utilizes a vacuum switch to detect the negative pressure at the beginning of each inspiration (approximately .1 - .25 cmH₂O due to inhalation). That, in turn, drives the solenoid for a time interval that corresponds to the flow rate selected on the rotary selector. At higher flow rates, the valve is open longer resulting in increased pulse volumes.

Atmospheric pressure compensation occurs automatically because the pressure side of the vacuum switch is open to atmosphere.

All PD1000A units are set on 2 liters per minute continuous flow from the factory. Changing to a different continuous flow rate involves changing the cannula fitting as outlined in Section VI, C.

The PD1000A delivers 16.5 cc O₂ per setting number (i.e. setting 2 = 33cc O₂).

INSTALLATION AND OPERATION

Installation

The PD1000A utilizes a rotary selector that has three modes for use by the patient. Those modes are: OFF, "PulseDose", and Continuous Flow. The PD1000A is battery operated and must be turned off by turning the Rotary Selector to the "OFF" position. The unit requires (2) "AA" batteries to operate in PulseDose mode (Figure 1). The battery does not become discharged as a result of not turning the Rotary Selector to the "OFF" position.

The PD1000A series conserving device is mounted on different size cylinders (Figure 2). It can be used on C, D, E, ML-6, M4 and M-6 size tanks at pressures between 500-2250 PSIG.

The rotary selector can be set to the prescribed flow rate. Verify that a pulse is being delivered at the leading edge of each inhalation. As the flow rate is increased, so is the duration of the pulse (Figure 3).

To use the PD1000A in the "Continuous Flow" mode, turn the Rotary Selector to the "CF" position (Figure 4). CF mode is not powered by the battery and can be used regardless of battery level. The "CF" flow rate is set at 2 liters per minute. This flow rate can be changed by changing cannula fittings (Figure 5). See Service Instructions Section VI, C.

Operation

When using the PD1000A in the PulseDose mode, the patient must breathe through the nose only. A standard nasal cannula must be used. Do not use pediatric or low-flow cannula. The cannula is connected to the cannula fitting. The cannula tubing can be up to 35 feet in length, but a 10 foot maximum is recommended to lessen the chances of the oxygen cylinder tipping over while in use.

Do not use on patients who can only mouth breathe. The PD1000A can only be used on patients capable of nose breathing and can not be used on mouth breathers.

Remaining battery life can be observed as the Pulse Indicator Lights illuminate with each breath (Figure 6).

If the Green Pulse Indicator Light illuminates with each breath, the batteries have sufficient power (8 or more hours).


If the Red Pulse Indicator Light illuminates with each breath, the batteries have between 4 – 8 hours of battery life remaining.

If the Red Indicator Light illuminates continuously, the batteries must be changed. The unit can be used on the continuous flow "CF" setting if no batteries are immediately available.

NOTE– Turn the Rotary Selector to the OFF setting prior to changing batteries.

NOTE– The oxygen cylinder will not last as long in continuous flow mode as it would in PulseDose mode.

The contents gauge indicates the approximate amount of oxygen left in the tank. It has a red area to emphasize when the tank needs to be changed (Figure 7).

MAINTENANCE PROCEDURES

Preventive Maintenance

NOTE– These tests only need to be performed after a repair or if the unit is going out on a new patient.


1. Install known good batteries.
2. Fill iFill cylinder to at least 500 PSI.
3. Verify that there are no leaks around the PD1000A and the oxygen tank.
4. Select any flow rate and simulate an inhalation through the nasal cannula while verifying that a "pulse" or bolus of oxygen is delivered with each simulated breath.
5. Position the Rotary Selector to the "CF" position and verify that the Continuous Flow Mode of operation is functioning.
6. Verify that the top and bottom covers are not cracked or broken.
7. Verify that the label can be read and is not damaged.
8. Turn the Rotary Selector to the "OFF" position.

Cleaning

The unit should be kept clean and free from moisture and dust. Clean the unit at least weekly by wiping with a dry, lint-free cloth. Avoid getting fluids or debris such as sand or dirt inside the oxygen connections. Do not immerse in water. Do not clean with a solvent based cleaning solution. Avoid dropping the unit or placing it in a position where it could topple or fall since this can damage the device. Whenever possible, use a padded carrying bag such as those listed in the Accessories section to carry the unit. This will help to protect it in case of a fall. The unit should be protected from extreme temperatures. Do not attempt any other maintenance.

Cylinder Recertification Testing

Cylinders will have two or three lines of markings stamped into the crown including the recertification test date. See examples below.



All US DOT marked iFill cylinders must undergo requalification testing every 5 years. Cylinders that are pi “π” marked, which is the official stamp mark required by the Transportable Pressure Equipment Directive (TPED) certifying that the gas cylinder conforms to RID/ADR for use throughout the European Union, will require recertification within 10 years as defined per Annex B of DIN EN 1802:2002 .

Recertification testing of cylinders must be done by a certified testing facility.

TROUBLESHOOTING

Please select from the types of problems listed below and follow the recommended action to logically troubleshoot and diagnose the problem.

PROBLEM	ACTION
Type I– The unit does not deliver a pulse with each inhalation, and the cylinder is pressurized.	1. If it is an international cylinder, check that the oxygen cylinder valve is open. 2. Verify that the batteries are good and of specified type. 3. If batteries are good, verify proper cannula connection and that patient is nose breathing. 4. If yes, remove top and bottom covers (Section VI, B). Use a digital volt meter and check the resistance across the Vacuum Switch (Figure 8) with Rotary Selector “OFF” when breath is simulated. 5. If the resistance doesn’t go to “0” when breath is simulated, change the vacuum switch. 6. If the resistance goes to “0” check the solenoid voltage (Figure 9) when breath is simulated and set to 6 LPM. 7. If no voltage is present, change the PC Board. 8. If voltage is present, change the solenoid.
Type II– The gauge is broken or not reading accurately.	Refer to Section VI to remove and replace the gauge.
Type III– The cannula fitting is broken.	Refer to Section VI, C to change the cannula fitting.
Type IV– Broken or missing battery door.	Replace battery door. (PD1000D-608)
Type V– Broken, damaged, or non-functioning rotary selector.	Remove covers and replace rotary selector of PC Board. Refer to Section VI, B.
Type VI– Unit won’t pulse.	Remove covers (Section VI, B). Test vacuum switch, PC Board, and Solenoid (Section VI, D).

SERVICE INSTRUCTIONS

WARNING

Ensure the iFill cylinder is empty before servicing the PD1000A. To empty a domestic cylinder, turn the rotary selector to the "CF" position and allow the oxygen to evacuate (Figure 10).

To empty a white, international cylinder, open the cylinder valve, turn the rotary selector to the "CF" position and allow the oxygen to evacuate (Figure 11).

Do not service the PD1000A until the cylinder is completely empty and the flow has stopped.

NOTE– Do not remove the PD1000A from the iFill cylinder. Most of the servicing and repairs can be performed while the device is attached to the cylinder. If the regulator or cylinder is damaged or requires removal for servicing, contact DeVilbiss Healthcare for return instructions.

Battery Removal

NOTE– Turn the rotary selector to the "OFF" setting and wait approximately 15 seconds prior to changing batteries.

Open the battery door by pushing back on the latch and lifting (Figure 12). Remove batteries ("AA") for replacement. Note the polarity when removing. The polarity is also indicated in the battery compartment.

Cover Removal/Gauge

Position the PD1000A face down so that the back cover is facing up. Remove the 5 cover screws (Figure 13) and remove the rear cover. Stand the unit on its side and (Figure 14) remove the circuit board/manifold assembly. Slide the front cover off the regulator. The Rotary Selector will also lift out when the front cover is removed (Figure 15). The gauge can then be removed from the regulator body. See Section VIII Parts and Illustrations.

Cannula Fitting Removal/Replacement

NOTE– When removing or installing the cannula fitting, orient the manifold so the fitting faces down to prevent any debris from falling into the manifold which could obstruct the continuous flow orifice.









Use a 5/16" open end wrench and remove the cannula fitting by turning it counter-clockwise. Replace with a new cannula fitting and torque to 10 inch-lbs.

Solenoid/Vacuum Switch Testing and Removal

Connect a piece of tubing to the cannula fitting so that a negative pressure can be created. Using a digital voltmeter, check the resistance across the vacuum switch as a negative pressure is created (Figure 8). If the resistance doesn't go to "0" when a breath is drawn, replace the vacuum switch.

If it goes to "0" but the solenoid doesn't open, check the solenoid voltage (3VDC). If there is no voltage, change the PC Board. If there is voltage, replace the solenoid (Figure 16).


To replace the vacuum switch, pull it apart from the manifold/solenoid and the PC Board terminals (Figure 17). The vacuum switch can then be replaced by pushing the manifold, vacuum switch, and PC Board together. No tools are required.

- | | | |
|----------------------------|--------------------|--------------------------|
| 1. Battery Indicator Light | 7. Cover Screws | 13. PC Board |
| 2. CF Position | 8. Slot for Gauge | 14. Front Cover |
| 3. Gauge | 9. Cannula Fitting | 15. Knob Opening |
| 4. Shut-Off Valve | 10. Solenoid | 16. Rotary Selector Knob |
| 5. Battery Door | 11. Manifold | 17. Solenoid Wires |
| 6. Battery Door Latch | 12. Vacuum Switch | 18. Solenoid Terminals |


PARTS AND ILLUSTRATIONS

Internal Parts


Part	Part Number	Qty
1. PC Board	PD1000D-607	1
2. Solenoid/Manifold	PD1000D-604	1
3. Vacuum Switch	PD1000D-605	1
4. Gauge	PD1000A-620	1
5. Regulator Body	PD1000D-606	1
6. Regulator Hose	PD1000D-620	1

External Parts

Part	Part Number	Qty
1. Cannula Fitting:		1
Aluminum- 2 lpm	PD1000D-612	
Gold- 3 lpm	PD1000D-613	
Green- 4 lpm	PD1000D-614	
Blue- 5 lpm	PD1000D-615	
Red- 6 lpm	PD1000D-616	
2. Cover Screws 6x.25	PD1000D-611	3
3. Cover Screws 6x1	PD1000D-610	2
4. Battery Door	PD1000D-608	1
5. Top/Bottom Cover	PD1000A-601	1
6. Rotary Selector	PD1000D-602	1
7. Gauge Protector	PD1000A-621	1
8. Fill Connector		

PNEUMATIC DIAGRAM

UNIT SPECIFICATIONS

Weight	26.3 ounces without battery 28 ounces with battery (PD1000A only)
Dimensions	4.13"L x 3.4"W x 4.75"H (10.49 cm L x 8.64cm W x 12.07 cm H) without cylinder
Power Supply	(2) Standard "AA" alkaline or NiMH.
NOTE – Batteries other than alkaline or NiMH are not recommended due to the capacity needed for operation and battery life of the unit.	
Operating Temperature Range.....	-10° to 40° C (14° to 104°F)
Operating Pressure Range.....	500 to 2250 PSIG (34 to 155 Bar) tank pressure
Operating Atmospheric Conditions.....	500 to 1020 millibar
Operating Humidity Range	0 to 95% R.H., non-condensing
Storage and Transportation Temperature Range.....	-40° to 70° C (-40° to 158° F)
Storage and Transportation Humidity Range	Up to 95% R.H., Non-condensing
Degree of Protection Against Ingress of Liquids	None
Degree of Protection Against Electric Shock	TYPE BF applied part
Modes of Operation.....	Continuous/Pulsed
Approval Body And Safety Standards:	Approved by CSA to: IEC 601-1;CAN/CSA-C22.2 No. 601.1-M90 as ordinary equipment
US Patents	5,755,224
Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.	

PARTS RETURN AND ORDERING POLICY

ALL DEFECTIVE COMPONENTS THAT ARE STILL UNDER WARRANTY MUST BE RETURNED TO THE FACTORY IN SOMERSET, PA WITHIN 30 DAYS AFTER SHIPMENT OF THE NEW COMPONENTS. IF THE COMPONENTS ARE NOT RECEIVED WITHIN THIS PERIOD, AN INVOICE WILL BE ISSUED TO YOUR ACCOUNT.

Before returning parts or units to the factory, call the DeVilbiss Healthcare Customer Service Department at 800-338-1988 or 814-443-4881 to obtain a return authorization number. Include in the package a note indicating the return authorization number along with your company name, address, phone number, and account number. The return authorization number should also be written on the outside of the package. To expedite your order for warranty or non-warranty parts, the following information should be given to the representative:

- Catalog Number
- Account Number
- Unit Serial Number
- Company name, address, and phone number

Ordering Information

When ordering components, instruction guides, or service manuals the following must be provided:

- Unit Catalog Number
- Part Number
- Unit Serial Number
- Quantity Required

Orders may be placed by calling:

Customer Service: 800-338-1988 • International Department: 814-443-4881 • Technical Service: 800-338-1988

PROVIDER'S NOTES

No routine calibration or service is required provided the device is used in accordance with the manufacturer's directions. Between patients wipe with a damp cloth having a maximum 5.25% Sodium Hypochlorite (Bleach) or 3% Hydrogen peroxide solution. Avoid getting fluids or debris such as sand or dirt inside the oxygen connections. Do not immerse in water.

DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	This device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker emissions	N/A	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	Complies	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5GHz	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz	N/A	
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential ±2kV common	N/A	
Power frequency magnetic field IEC 61000-4-8	3 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interrupts and voltage variations on power supply input lines IEC 61000-4-11	>95% dip 0.5 cycle 60% dip 5 cycles 70% dip 25 cycles 95% dip 5 secs.	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

This device has been tested to and meets the EMC requirements of EN60601-1-2. Do not place the device near other equipment or devices that create or attract electromagnetic fields. Examples of such equipment are defibrillators, diathermy equipment, CB radios, microwave ovens, etc. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

DeVilbiss Healthcare LLC
100 DeVilbiss Drive
Somerset, PA 15501-2125
USA
800-338-1988 • 814-443-4881

DeVilbiss Healthcare Ltd
Unit 3, Bloomfield Park
Bloomfield Road
Tipton, West Midlands DY4 9AP
UNITED KINGDOM
+44 (0) 121 521 3140

DeVilbiss Healthcare Pty. Limited
15 Carrington Road, Unit 8
Castle Hill NSW 2154
AUSTRALIA
+61-2-9899-3144

DeVilbiss Healthcare SAS
13/17, Rue Joseph Priestley
37100 Tours
FRANCE
+33 (0) 2 47 42 99 42

DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim
GERMANY
+49 (0) 621-178-98-0

DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA
800-338-1988 • 814-443-4881 • www.DeVilbissHealthcare.com

DeVilbiss® and PulseDose® are registered trademarks of DeVilbiss Healthcare.

© 2015 DeVilbiss Healthcare LLC. 03.15 All Rights Reserved.

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