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Table of contents

2. Definition of the target groups.	32	7. Operation.	45
2.1 Users	32	7.1 Function test.	45
2.2 Servicing personnel	32	7.2 Beginning operation	45
3. For your safety and the safety of your patients	32	7.3 Ending operation	46
3.1 General safety information	32	8. Faults, trouble-shooting	46
3.2 Carefully observe the instruction for use	32	8.1 General.	46
3.3 Servicing	32	8.2 Replacing the oscillating quartz	47
3.4 Accessories	33	8.3 Replacing fuses	48
3.5 Product-specific safety information	33	9. Treatment.	48
3.6 Patient safety	33	9.1 Disassembly	48
3.7 Duty to report serious incidents	33	9.2 Execution instructions	48
3.8 Information regarding electromagnetic compatibility	33	9.3 Performing cleaning.	49
3.9 Checking the delivery.	33	9.4 Manual cleaning of the quartz	50
3.10 Symbols	34	9.5 Cleaning run-on kits.	50
4. Safety information.	34	10. Checks / repairs	51
5. Device description.	36	10.1 General.	51
5.1 Device description.	36	10.2 Checks	51
5.2 Intended use / indication.	36	10.3 Repairs	51
5.3 Application area and conditions of use	36	11. Delivery program.	52
5.4 Contraindications	37	11.1 Standard delivery (2010-41-0-000, Class IIa)	52
5.5 Reasonably foreseeable misuse	37	11.2 Optional accessories	52
5.6 Device components.	37	11.3 Spare parts and consumables	52
5.7 Function description / frequently used functions	38	11.4 Interface description	53
6. Commissioning	40	12. Disposal	53
6.1 Preparation	40	13. Technical data	54
6.2 Establishing the power supply	40	13.1 Environmental conditions	54
6.3 Inspection	41	13.2 Performance data	54
6.4 Preparing the USV 2010	41	14. Electromagnetic compatibility	55
6.4.1 Nebulisation directly from the nebuliser chamber.	41	14.1 Electromagnetic emissions	55
6.4.2 Nebulisation with commercially available disposable sterile water systems	42	14.2 Electromagnetic interference on other equipment	55
6.4.3 Nebulisation with run-on system (optional)	43	14.3 Electromagnetic interference immunity.	56
6.5 Fitting accessories.	43	14.3.1 Guidelines and manufacturer's declarations - electromagnetic interference immunity.	56
6.5.1 Installing the table-top model (optional accessory).	43	14.3.2 Guidelines and manufacturer's declarations - electromagnetic interference immunity, non-life-supporting equipment.	57
6.5.2 Installing the wall-mounted model accessory kit (optional accessory)	44	14.4 Recommended safety distances	57
6.5.3 Installing the carriage accessory kit (optional accessory)	44		
6.5.4 Transport with the carriage.	45		

1. Definition of the safety information

WARNING!

Information with the word **WARNING** warns of a dangerous situation that can lead to death or serious injuries.

CAUTION!

Information with the word **CAUTION** warns of a situation that can lead to minor or moderate injuries.

ATTENTION!

Information with the word **ATTENTION** warns of a situation that can lead to damage to property or the environment.

NOTE

Provides additional information.

2. Definition of the target groups

The product may be used, installed, treated or serviced exclusively by users from the defined user groups, as well as maintenance personnel. These target groups must possess the training and specialist knowledge required to use, install, set up or service the product.

2.1 Users

Users are persons who use the product in accordance with the intended use/indication.

2.2 Servicing personnel

Servicing personnel are persons who are responsible for servicing the product.

Servicing personnel are persons who have been trained on servicing medical devices and who install, treat or service the product.

3. For your safety and the safety of your patients

3.1 General safety information

The following text marked with **WARNING** or **ATTENTION** is general safety information for the operation of the medical device. Text marked with **WARNING** or **ATTENTION** pertaining to special parts or functions of the medical device appears in the affected sections of this instruction for use or in the instruction for use of another product that is used together with this medical device.
Keep the user manual easily accessible near the device!

3.2 Carefully observe the instruction for use

WARNING!

Each application of the medical device requires precise knowledge and observation of all chapters of this instruction for use. The medical device must be used exclusively for the purpose specified under intended use. All safety information marked with **WARNING** or **ATTENTION** in this instruction for use and the information on the medical device labels must be carefully observed. A failure to observe this information is considered misuse of the medical device contrary to the intended use.

3.3 Servicing

ATTENTION!

Danger of a malfunction of the medical device and patient endangerment.

The medical device must be subjected to regular inspection and maintenance by maintenance personnel. Danger of a malfunction of the medical device and patient endangerment with a failure to observe the above. Observe information regarding maintenance / repair.

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3.4 Accessories

The manufacturer hereby declares that they are only responsible for the safety, reliability and serviceability of the device if original spare parts are used and if device repairs are carried out exclusively by persons who have been authorised by the manufacturer to do so. The statutory warranty period applies. In the event of warranty claims, please contact the authorised dealer.

⚠ WARNING!

Only the accessories listed in the order list have been tested and approved for use with the medical device. It is therefore strongly recommended that the medical device be used in conjunction with these accessories only. Otherwise the correct function of the medical device may be at risk and no guarantee is provided.

3.5 Product-specific safety information

⚠ WARNING!

Danger of malfunction!

Unauthorised modifications of the medical device result in malfunctions and may lead to dangerous situations.

It is prohibited to modify this medical device without the manufacturer's permission.

3.6 Patient safety

With the equipment of the medical device, the associated documentation and the medical device labels it is assumed that the medical device is acquired and used exclusively by specialists and that the basic functions of the medical device are known to the user. The instructions for use as well as the safety information marked **WARNING** and **ATTENTION** are therefore essentially limited to the special features of the medical device. This instruction for use does not contain any information regarding hazards that are obvious to professionals, the consequences of incorrect use of the medical device and possible negative effects on patients with various underlying diseases.

3.7 Duty to report serious incidents

If serious incidents should occur in conjunction with this device, the incidents must be reported to the following bodies:

- The manufacturer.
- The authority responsible in the country in which the incident took place.

3.8 Information regarding electromagnetic compatibility

General information regarding electromagnetic compatibility (EMC) in accordance with international EMC standard IEC 60601-1-2.

Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC) and must be installed and commissioned according to the EMC information provided.

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Portable and mobile high frequency communication devices can affect medical electrical devices.

3.9 Checking the delivery


When the device is delivered, check that all components are present and in faultless condition. Take the device out of the packaging and first check that the delivery matches with the order and that it is complete.

⇒ „11.1 Standard delivery (2010-41-0-000, Class IIa)“ (Seite 52).

In the event of complaints, please contact the authorised dealer. Devices are not delivered sterile.

3.10 Symbols

The following symbols are on the device:

	Observe the instruction for use
>ABS<	Material designation of the plastic
	CE-mark and number of the notified body according to Medical Device Directive 93/42/EEC
	Standby push-button
	Mist volume
	Air volume
	Timer
	Device with application part, type B
	Law restricts this to sale by or on the order of a physician
	Quantity
	Article number
	Serial number
	Revision index
	Date of manufacture (year-month-day)
	Manufacturer
	Protect from moisture
	AC voltage
IP 22	Dripping water protected
	Medical device
	Device fuse
	Labelling in accordance with EC Directive 2002/96/EC (Waste Electrical and Electronic Equipment Directive) Symbol for "Do not dispose of product at municipal collection points for used electrical equipment"

4. Safety information

WARNING!

The ventilation slots must not be covered by objects, otherwise the device will not operate correctly and the aerosol nebulisation will be destroyed.

WARNING!

Risk of infection! Observe hygiene regulations!

- The applicable hygiene rules must be complied with to avoid infection or bacterial contamination. Observe the intended use of the bacteria filter. Always wear gloves during work.
- Contaminated device parts may endanger the patient's health. Before use, the device must be treated in accordance with the operational hygiene regulations.

WARNING!

Danger of injury!

The connection for the optional heating tube on the side of the device must not be touched at the same time as the patient.

ATTENTION!

Damage!

The device must be checked for damage (housing, etc.) and safety-related defects every time before use. Damaged devices must be serviced by an authorised party before reuse. All work that requires the use of tools must be performed by the manufacturer or a party authorised by the manufacturer.

ATTENTION!

Power supply!

The mains power supply must concur with the specifications on the type plate. Unplug the mains plug before cleaning. Prevent any liquid penetrating the device interior.

Safety information

⚠ ATTENTION!

Liability!

The manufacturer accepts no liability and provides no warranty for other uses that deviate from the instruction for use, or for the unauthorised use of non-original spare parts. Repairs may only be carried out by the manufacturer or authorised specialists using original parts and in compliance with the safety precautions applicable to medical electrical devices.

⚠ ATTENTION!

Function testing!

The user must check the functionality and suitability of the application parts for the respective purpose as well as their biocompatibility.

5. Device description

5.1 Device description

The device consists of a washable plastic housing with internal electronic components and an attachable nebuliser chamber. The device is started up and controlled with push-buttons. The base of the nebuliser chamber contains an oscillating quartz, which causes the liquid in the container to nebulise through its high-frequency oscillations. A blower located in the housing conveys the aerosol through the patient tube to the desired breathing area of the patient. The aspirated air is passed through a bacteria filter before nebulisation. The device is designed for a service life of 10 years.

Three different set-up variants are provided for different installation options of the ultrasonic nebuliser:

- Table-top version:
The device stands on a table or on a horizontal and fixed support surface
- Wall-mounted version:
The device is fastened to a device rail 25x10mm (DIN)
- Carriage version:
The device is fastened to the stand tube of a carriage

5.2 Intended use / indication

The ultrasonic nebulisers are exclusively intended for humidifying the respiratory air of spontaneously breathing patients from the age of 16 years and for aerosol therapy.

Aerosol therapy should only be carried out according to the instructions of a doctor.

Inhalation for humidification of the respiratory tract using distilled water or brine solution 0.9%.

With this application the aerosol causes a mucus solution in the respiratory tract and thus promotes the expectoration of the secretion. In addition, the dilution of the bronchial secretion supports the body's own cleaning function of the lungs.

All other uses are deemed non-intended use. The manufacturer shall not be liable for damage caused by non-intended use.

5.3 Application area and conditions of use

The medical device shall be used exclusively in rooms designated for clinical practice.

WARNING!

Danger of explosion!

- The medical device is not approved for operation in areas in which oxygen concentrations of more than 25% by volume, flammable or explosive gas mixtures may occur.
- Do not use the device together with flammable or explosive liquids.

WARNING!

Patient endangerment!

- Do not use the ultrasonic nebuliser on patients requiring ventilation.
- Do not install ultrasonic nebulisers in breathing systems.
- Do not use the ultrasonic nebuliser with a respiratory mask, tracheal cannula.
This can lead to fluid accumulations in the respiratory organs.

WARNING!

Risk of infection!

Do not use ultrasonic nebulisers without a bacteria filter.

The field of application of the ultrasonic nebuliser USV 2010 is a clinical environment, medical practice, senior citizens' or nursing home and in the military medical service by qualified personnel trained in the specific application. The ultrasonic nebuliser is intended for short-term use.

Due to the alveolar particle size of the mist, special care must be taken with regard to hygiene when preparing and performing inhalation.

The air flow is only sufficiently filtered if a suitable bacteria filter is used.

Replace the bacteria filter after 24 hours of use or according to the special instructions of the respective filter manufacturer.

Clean reusable tubing daily or when changing patients.

Replace disposable tubes and mouthpieces daily or when changing patients.

Device description

5.4 Contraindications

- Coughing, bronchospasm or allergic reactions.
- Inability of the patient to cooperate accordingly.

5.5 Reasonably foreseeable misuse

⚠ WARNING!

Observe application area!


The device must not:

- be covered, thrown, pushed open or dropped,
 - be used in the ventilation of patients,
 - be operated without bacteria filter or liquids,
 - be operated by children when unsupervised.
-
- Do not operate with flammable or explosive liquids or in hazardous areas.
 - Do not insert objects in the housing.


5.6 Device components

- 1 Connection for bacteria filter
- 2 Ultrasonic nebuliser
- 3 Tube approx. 300 mm
- 4 Bacteria filter
- 5 Tube approx. 1,200 mm
- 6 Nebuliser chamber
- 7 Nebuliser chamber cover
- 8 Run-on kit
- 9 Quartz module
- 10 Mains cable
- 11 Nebuliser chamber, complete

- 12 Main switch
- 13 Mounting bracket
- 14 Device socket
- 15 Voltage selector
- 16 Connection socket for heating tube (accessories)


- 17 Retainer part for nebuliser chamber
- 18 Standby push-button
- 19 Mist volume push-button
- 20 Air volume push-button
- 21 Timer push-button

5.7 Function description / frequently used functions

In the ultrasonic nebuliser USV 2010, high frequencies are generated by an oscillating quartz which converts liquids inside the nebuliser chamber into ultra-fine mist. The mains cable (10) is plugged into the device socket (14). Pay attention to correct connection values in accordance with the information plate. The ultrasonic nebuliser is switched on with the main switch (12) located on the back of the unit. The aspirated air is guided through the bacteria filter (4) via the 300 mm tube into the nebuliser chamber (11). The mist of liquid produced by means of oscillating quartz inside the nebuliser chamber (11) is conducted with the generated air flow through the patient tube (5) to the vicinity of the patient's respiratory organs.

Standby push-button

- ➡ Switch on the device with the standby button [1] on the front of the device.
- Green LED [2] lights up/field designation: 
- The device starts nebulising.


NOTE:

In case of faults, the red LED [3] lights up.

Device description

Mist volume push-button

The push-button is used to regulate the amount of mist produced per time in the nebuliser chamber.


Increase mist volume

- ☞ Press the button (+) [1].
- One further yellow LED [2] lights up respectively.

Reduce mist volume

- ☞ Press the button (-) [1].
- One yellow LED [2] goes out.

Air volume push-button


Increase air volume

- ☞ Press the button (+) [1].
- One further yellow LED [2] lights up respectively.

Reduce air volume

- ☞ Press the button (-) [1].
- One yellow LED [2] goes out.

Setting the timer:

- ☞ Pressing the button (+) or (-) sets the desired time on the timer, see [2].
Time can be set up to 60 min.
- When the timer is active, the yellow control lamp lights up and shows the set time.
 - Once the set time has elapsed, the device switches off automatically:
 - The red lamp [3, see information regarding standby push-button] lights up and an acoustic warning tone sounds. After approx. 2 minutes, the warning tone ceases.
 - Switching on again can only take place after switching off with the standby push-button

NOTE:

When using sterile water capsules without a run-on system, do not exceed a max. setting time of 15 minutes.

6. Commissioning

6.1 Preparation

Place the device upright on a level, firm surface. The correct function of the overflow protection is only guaranteed if the device is in vertical position. If placed on a soft surface (e.g. pillow, mattress) there is a risk that the ventilation slots will be covered and the device will overheat.

6.2 Establishing the power supply

⚠ WARNING!

Danger of electric shock!


Connecting devices to additional sockets can lead to an increase in leakage current. If the protective conductor of one of these devices fails, the leakage current may exceed the permissible values. Only establish connections with the approval of the respective device manufacturer. Have servicing personnel check the leakage current. If the permissible value is exceeded, do not use the additional sockets, but separate mains sockets on the wall.

⚠ ATTENTION!

Device out of order!

Before initial commissioning, the device voltage must be set to the available mains voltage using the voltage switch on the back of the device. 230V is pre-set.


Selecting mains voltage

- ✎ Carefully unlock and remove the fuse link (A) at the device socket with a screwdriver.
- ✎ Remove the voltage switch (B) from the fuse link and select the desired mains voltage by turning it. Replace the fuse link by pressing it until it audibly clicks back into place.
- The required mains voltage (110 V/ 127 V/ 230 V) is now visible in the window (C).

NOTE:

Observe type plate on the underside of the device.

- ✎ Insert the plug of the mains cable (D) in the device socket and in a mains socket with earthing contact.

Commissioning

6.3 Inspection

Only parts of the device that are in perfect condition guarantee the functionality of the ultrasonic nebuliser. Therefore, the device parts must be thoroughly checked and inspected every time before assembly.


Inspect and check:

- whether the mains cable is free of damage,
- whether the plastic and rubber device parts (e.g. tube, nebuliser chamber) are in perfect condition and that no damage due to ageing is apparent,
- whether the device has been correctly cleaned and is clear of residues or contaminants.

6.4 Preparing the USV 2010


6.4.1 Nebulisation directly from the nebuliser chamber

Inserting the nebuliser chamber

- ☞ Fill the nebuliser chamber (2) with nebuliser liquid up to the max. mark [2]. Do not fill beyond the max mark [2].
- ☞ Insert the nebuliser chamber (2) in the holder of the USV 2010 [1]. Observe preferred position.
- ☞ Press cover onto the nebuliser chamber.


Fitting the bacteria filter

- ☞ Plug the bacteria filter [1] onto the blower connection [2] on the housing.

Pay attention to the flow direction according to the manufacturer's specifications.


Fit 300 mm tube

A 300 mm tube (1) connects the bacteria filter (2) and the cover of the nebuliser chamber (3).

- ☞ Place the short tube [1] on the nozzle of the bacteria filter [2] and the nozzle of the lid [3].

Fit 1200 mm tube

- ☞ Place the long tube [1] on the second nozzle of the cover [2] and press it into the clamps of the articulated arm. Guide the tube upwards so that condensate flows back into the nebuliser chamber.

Secure the cover of the nebuliser chamber with the retainer.

- ☞ Prepare the sterile water system in accordance with the corresponding instruction for use and insert it in the nebuliser chamber in place of the cover.
- ☞ Connect the tubes and secure the sterile water system with the retainer.
- ☞ Change the contact water after each use, although at least 1x per day.

NOTE:

Sterile water capsules are available e.g. in pharmacies or in medical supply stores.

Ensure that the nebuliser chamber and the sterile water system have a sufficient liquid level.


If this is not the case, the sterile water capsule may be damaged (loss of sterility) or the nebuliser chamber may be damaged. Always secure sterile water capsules with the retainer (25).

⚠ WARNING!

Danger of injury!

When the sterile water capsule is empty, the water in the nebuliser chamber can heat up to over 70°C.

6.4.2 Nebulisation with commercially available disposable sterile water systems


- ☞ Fill the nebuliser chamber (2) with distilled water until the sterile water system [2] (e.g. Hudson RCI) is wetted up to mark for the contact water (hatched area [1]). The distilled water in the nebuliser chamber serves to transfer the vibrations to the sterile water system.

Commissioning

6.4.3 Nebulisation with run-on system (optional)

Run-on systems are used to provide a permanent supply of liquids to fill the nebuliser chamber.

- ✎ Fit the prepared run-on system [3] in the nebuliser chamber in place of the cover.
- ✎ Connect the tubes and secure the run-on system with the retainer.
The height difference between the nebuliser chamber and container (bag/bottle) must not exceed 1000 mm.
- ✎ Tap the container with a transfer kit, e. g. Intrafix Air, brand: B.Braun. Do not bend the cable.
Ensure sterile venting of the container.

Function test of the run-on system

After connecting the container, check the function of the float:

- Fluid runs into the nebuliser chamber.
- Intake stops automatically when the liquid level reaches its level (range between min. and max. marking).
- ✎ If a fault occurs, clean the float valve and check its function again. Then put the ultrasonic nebuliser into operation.

Ending operation of the run-on system

After operation, lock the intake of the container, e. g. by closing the knurled screw on the transfer kit.


6.5 Fitting accessories

⚠ ATTENTION!

Installation on the base unit must be carried out according to the operating instructions of the base unit on which this medical device is used. Check the safe connection to the base unit system.

- ✎ Carefully follow the instructions for use and installation.

6.5.1 Installing the table-top model (optional accessory)

- ✎ Insert the stand tube (300 mm long) into the rear retainer and fix it with the locking screw.
- ✎ Slide down retainer (to secure the cover of the nebuliser chamber) and articulated arm onto the stand tube and adjust.

Avoid a dip in the patient tube.


⇒ „6.4 Preparing the USV 2010“ (Seite 41).

NOTE:


- ✎ Now place the adjusting ring of the bottle holder (24) on the stainless steel stand tube and secure it in the desired position with the locking screw.
- ✎ Hang a container with the liquid to be nebulised on the bottle holder.

The air inlet on the container must be opened so that the infusion solution flows sufficiently into the nebuliser chamber.

6.5.2 Installing the wall-mounted model accessory kit (optional accessory)


- 22 Stainless steel stand tube 25 x 600 mm
- 23 Articulated arm with adjusting ring
- 24 Bottle holder with adjusting ring
- 25 Retainer with adjusting ring
- 26 Holding claw for device rail


- ☞ Attach the holder and stand tube (600 mm long) to a device rail (25x10 mm) on the wall.
- ☞ Fit the USV 2010 with the rear retainer on the stand tube and fasten it with the locking screw.
- ☞ Slide down retainer (to secure the cover of the nebuliser chamber) and articulated arm onto the stand tube and adjust.

6.5.3 Installing the carriage accessory kit (optional accessory)


- 27 Stainless steel stand tube
- 28 Articulated arm with adjusting ring
- 29 Bottle holder with adjusting ring
- 30 Retainer with adjusting ring
- 31 3-foot carriage


- ☞ Push the locking ring (A) onto the stand tube and screw it tight.
- ☞ Fit the USV 2010 with the rear retainer on the stand tube of the carriage and fasten it with the locking screw.
- ☞ Slide down retainer (to secure the cover of the nebuliser chamber), articulated arm, bottle holder (for commercially available top-up systems) onto the stand tube and adjust.

Operation

6.5.4 Transport with the carriage

☞ When locking the brake wheels [B], make sure that the USV is aligned between the brake wheels.

⚠ WARNING!

Risk of toppling!

For transport within the clinic or hospital, push the USV down to the pre-assembled stop (max. 70 cm above the floor) and fold it into the jointed arm. Guide the device during transport.

⚠ ATTENTION!

Risk of crushing during assembly, commissioning and use of the articulated arm.

7. Operation

7.1 Function test

⚠ WARNING!

The ultrasonic nebuliser is used for the treatment of patients exclusively according to its intended use. A performance limitation of the device can lead to complications during treatment.

Before each use, the ultrasonic nebuliser must be subjected to a complete functional test. Make sure that

- the bacteria filter is installed and functional,
- the tubes and nebuliser chamber are assembled and functional and there are no residues in these,
- the tubes and the nebuliser chamber are secure, leak-tight and no mechanical forces are applied to the tubes; the tubes must not be kinked,
- device parts made of plastic or rubber are in perfect condition and show no signs of ageing,
- the mains cable is free of damage,
- the nebuliser chamber is sufficiently filled,
- the retainer rests on the respective cover.

7.2 Beginning operation

Only carry out aerosol therapy according to the instructions of a doctor!

- ☞ Aim the open end of the long tube toward the mouth and nose area.
- ☞ Maintain a minimum distance of 20 cm from the patient.
- ☞ Set the toggle switch on the rear to "I" (on).
- ☞ Set the desired nebuliser output.
- ☞ Set the desired air flow.
- ☞ Ensure an adequate liquid level.


Do not allow the nebuliser chamber to run dry!

If the device runs dry, it switches off automatically and a signal tone sounds. After approx. 2 minutes, the signal tone ceases.

If the red lamp indicates a malfunction:

- ☞ Switch off the device with the standby push-button.
- ☞ Switch (off) the main switch on the rear to "0".
The red lamp goes out.
For troubleshooting:
⇒ see „8. Faults, trouble-shooting“ (Seite 46).
- ☞ Switch the main switch back to "I" (on) and switch on the device with the standby button.

7.3 Ending operation

Press the standby push-button [1].

All LEDs are extinguished and the device stops all functions.

- ☞ Switch the main switch on the rear to "0" (off).
- ☞ Close the tube clamp of the top-up system.
- ☞ If not used for an extended period of time (more than 60 minutes) unplug at the mains plug.

8. Faults, trouble-shooting

8.1 General

The ultrasonic nebuliser was subjected to a thorough quality inspection at the factory. However, if a fault should occur you may be able to rectify this yourself if you observe the following instructions.

On the left side of the control panel, there is a fault indicator on the standby push-button. If a fault occurs or the set time (end of timer) is terminated, this lights up red.

Once the set time has come to an end:

- ⇒ „5.7 Function description / frequently used functions“ (Seite 38).

If a fault occurs, first switch off the device.

- ☞ To do this, press the main switch (2).
If possible, eliminate the cause of the fault and switch the device on again.

The faults are described as follows:

What fault is present?

Cause of the fault.

- ☞ Fault remedy.
⇒ Reference to corresponding chapter.

The device does not switch on / the red lamp lights up.

The plug is not correctly inserted in the device or socket.

- ☞ Check the plug is sitting correctly in the device or socket.

The main switch on the rear of the device is not set to on (I).

- ☞ Switch the main switch to on (I).

Mains fuse missing, defective or incorrect.

- ☞ Check the mains fuse.

- ☞ Replace the fuse.

⇒ „8.3 Replacing fuses“ (Seite 48).

The thermal switch has triggered.

- ☞ Leave the device to cool down. The device is functional again after approx. 10 minutes.

The electronics are defective.

- ☞ Have the device repaired.
-

Faults, trouble-shooting

The device does not work / the red lamp lights up

The nebuliser chamber is inserted incorrectly.

- ☞ Insert the nebuliser chamber correctly. Note the preferred position.

The nebuliser chamber is empty.

- ☞ Top up the nebuliser liquid.

The nebuliser chamber is leaking.

- ☞ Screw the oscillating quartz firmly into the nebuliser chamber.

The oscillating quartz is defective.

- ☞ Replace the oscillating quartz.

The thermal switch has triggered.

- ☞ Leave the device to cool down. The device is functional again after approx. 10 minutes.

The electronics are defective.

- ☞ Have the device repaired.

Misting too low despite setting to maximum nebuliser output.

The nebuliser liquid is off the scale.

- ☞ Adjust the level of the nebuliser liquid to the scale (top up or remove nebuliser liquid).

The resistance of the filter is too high.

- ☞ Replace the filter.

There is condensate in the tube.

- ☞ Drain the condensate.

Guide the tube such that condensate flows back into the nebuliser chamber.

The tubes or the cover of the nebuliser chamber are leaking.

- ☞ Put the cover on tightly and secure it with a retainer.
- ☞ Push the tubes firmly onto the grommets.

No nebulisation is taking place.

The oscillating quartz has failed.

- ☞ Change the quartz module.

Contact water missing.

- ☞ Top up the contact water.

- ☞ Check the chamber for leak-tightness.

An accumulation of condensation water in the nebuliser tube prevents aerosol flow.


- ☞ Drain the condensate water.

⇒ „6.4 Preparing the USV 2010“ (Seite 41).

The bacteria filter or tubes are incorrectly fitted.

⇒ „6.4 Preparing the USV 2010“ (Seite 41).

8.2 Replacing the oscillating quartz

- ☞ Screw a new oscillating quartz (C) into the nebuliser chamber until the oscillating quartz is firmly seated and the edge of the nebuliser chamber is flush with the lower edge of the quartz module.

If it is not possible to remedy a fault:


- ☞ Contact the device supplier.

8.3 Replacing fuses

⚠ WARNING!

Electric shock!

Disconnect the device from the mains before replacing the fuses.

- 1 Device socket
- 2 Fuses
- 3 Locking mechanisms
- 4 Fuse insert

⚠ WARNING!

Observe fuse type!

Only fuses of the following type may be used:
T1.0 A H 250V 5x20 mm

To replace the fuses in the device socket, proceed as follows:

1. Unplug the mains cable from the device socket (1).
2. Carefully unlock the side locking mechanisms (3) on the device socket (1) with a small screwdriver and remove the fuse insert (4).
3. Remove defective fuses (2), replace with new fuses.
4. Insert fuse insert (4) and press into place until it audibly latches.
5. Check secure seating of the fuse insert (4).

9. Treatment

9.1 Disassembly

- ✎ Unplug mains plug and device plug.
- ✎ Unplug the tube heating plug.
- ✎ Remove tubes and bacteria filter.
- ✎ Slide the retainer up and remove the cover from the nebuliser chamber.
- ✎ Pull the nebuliser chamber out of the connection and pour out liquid from the nebuliser chamber. Observe the hygiene regulations of the hospital!
- ✎ Dispose of disposable items (disposable tubes and covers, containers of sterile water systems) in an environmentally friendly manner. Observe hygiene regulations.
- ✎ Do not disinfect or sterilise filters.
- ✎ Replace and dispose of filters according to the manufacturer's instructions.

9.2 Execution instructions

The individual parts of the ultrasonic nebuliser can be cleaned and disinfected according to the specifications in the following table. When doing so, please observe the following points:

- Before disinfection the parts must always be cleaned and then dried.
- Dirt particles can become encapsulated and cause equipment parts to be non-sterile after disinfection. Therefore, all parts of the device must be thoroughly cleaned before disinfection and stubborn dirt must be removed with a cloth and detergent.
- When using disinfectants, please follow the instructions of the manufacturer and the hygiene specialist.

⚠ WARNING!

Danger of injury!

Disinfectant can contain substances that are hazardous to health and that can cause irritation upon contact with skin or eyes. Protect your skin and eyes and observe the hygiene rules when working with disinfectants.

Treatment

⚠ ATTENTION!

Material changes!

Almost all components of the ultrasonic nebuliser are made of plastics. Solvents, some disinfectants or cleaning agents can dissolve the plastics or cause stress cracks. Therefore do not use any agents containing hydrocarbons or alcohol to clean the surfaces.

⚠ ATTENTION!

Function testing!

After cleaning or disinfection, always check that all parts of the device are in good working order.

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NOTE:

Discoloration of the plastic parts may occur when using non-colourfast wrapping cloths.

NOTE:

After cleaning, disinfection and complete drying, always check the rubber seal of the quartz module carefully for possible damage.

9.3 Performing cleaning

General

Rinse the nebuliser chamber and nebuliser chamber cover with distilled water after each use. Disinfection must be carried out with standard medical disinfectants.

The quartz module remains firmly screwed into the nebuliser chamber.

Rinse all individual parts with distilled water after each use.

Disposable item

Intended for single use only.

Disposable items must not be re-treated, there is a risk of infection.

- Bacteria and virus filter:
Replace after 1 week at the latest or according to special instructions from the manufacturer.
When using a closed sterile water system, it is recommended to change the filter after 48 hours of operation or every 14 days.

Cleaning

In warm water with mild household cleaner:

- Nebuliser chamber (transparent)
- Nebuliser cover (transparent)

Wipe with a damp cloth:

- Housing
- Mains cable
- Accessory kits

Disinfection

After the prescribed exposure time (according to the manufacturer's instructions), components must be cleaned of disinfectant residues with a damp cloth and then dried.

Immerse in diluted solution:

- Nebuliser chamber (transparent)
- Nebuliser chamber cover (transparent)

Wipe disinfection:

- Housing
- Mains cable
- Accessory kits

Cleaning in a washing machine

Rinse cycle up to 95° C:

- Nebuliser chamber (transparent)
- Nebuliser chamber cover (transparent)

Not permissible:

- Housing
- Mains cable
- Accessory kits

Autoclaving

Attention, hot steam accelerates the natural ageing of plastics. The function of the device parts can be impaired by the change in material.

Up to **121° C**, max. 20 min

- Nebuliser chamber (transparent)
- Nebuliser chamber cover (transparent)

Not permissible:

- Housing
- Mains cable
- Accessory kits

Examples of surface disinfectants:

- Incidin® Plus *
- Incidin® Perfect *
- Antiseptica Kombi surface disinfectant *
- Antifect® FF *
- B 10 *

Examples of instrument disinfectants:

- Cidex® AF
- Gigasept® FF *

* registered trademarks of the respective manufacturers

The lifespan of the nebulizer chamber with lid and the quartz module is limited to 30 autoclaving cycles. The aforementioned example disinfectants are recommendations from the device manufacturer. The recommendation of the mentioned active ingredient bases (pH 4.5-8) applies. This does not release the user from the need to consult a hygiene specialist in individual cases or to test the components to be disinfected for compatibility with the active ingredients.

9.4 Manual cleaning of the quartz

- ☞ To clean, remove the nebuliser chamber cover (7) from the nebuliser chamber (6) and unscrew the quartz module (9) from the nebuliser chamber (6). Reassembly takes place in reverse order.

9.5 Cleaning run-on kits

- ☞ Please refer to the separate instruction for use for instructions.

10. Checks / repairs

10.1 General

CAUTION!

Danger of harm to health!

The ultrasonic nebuliser is used for the treatment of patients. The device or parts of it could be contaminated. Before returning the device for inspection or repair, the nebuliser chamber, the run-on set, the bacteria filter and all tubes must therefore be removed, and the device must be cleaned and disinfected.

10.2 Checks

Maintenance work by the operating personnel beyond the usual care measures is not necessary. To ensure the functionality of the product, visual and functional tests must be carried out before use. To ensure the operational safety of the product, the availability of all functions and an increased service life, we recommend an annual inspection by the manufacturer. The scope of the inspection depends on a device-specific performance plan. The safety test in this regard is carried out in accordance with VDE 0750-1. The distributor recommends an annual "STK" (technical safety check).

10.3 Repairs

Disclaimer

All liability is excluded in the event of interventions by unauthorised persons. Repairs must be carried out by the manufacturer exclusively.

If defects are detected, the product shall no longer be used. In this case, send the device to the manufacturer for inspection or repair.

As soon as the reusable accessories show excessive signs of use, the respective component must be replaced.

Never carry out inspection or repair work whilst a patient is being treated with the device.

11. Delivery program

11.1 Standard delivery (2010-41-0-000, Class IIa)

comprising:

- Ultrasonic nebuliser base unit
- Mains cable
- Sterilisable nebuliser chamber complete with nebuliser chamber cover and quartz module
- Disposable corrugated tube approx. 300 mm
- Disposable corrugated tube approx. 1150 mm
- Bacteria filter

11.2 Optional accessories

The following accessories are not included in the delivery and must be ordered separately:

Carriage accessory kit (2002-00-5-900, Class I)

comprising:

- 5-foot carriage, chrome-plated, 50 mm castors,
- stainless steel stand tube, articulated arm with adjusting ring,
- bottle holder with adjusting ring (stainless steel),
- retainer with adjusting ring (stainless steel)

Table-top model accessory kit (2002-01-3-000, Class I)

comprising:

- stainless steel stand tube 25 x 600 mm, articulated arm with
- adjusting ring, bottle holder with adjusting ring (stainless steel),
- retainer with adjusting ring (stainless steel)

Wall-mounted model accessory kit (2002-01-4-000, Class I)

comprising:

- stainless steel stand tube 25 x 600 mm, articulated arm with
- adjusting ring, bottle holder with adjusting ring (stainless steel),
- retainer with adjusting ring (stainless steel), holding claw for device rail 25 x 10 mm (DIN)

11.3 Spare parts and consumables

2002-02-6-900, Class IIa

Supply tube for run-on kit

For ultrasonic – supply tube for run-on kit

2002-02-1-955, Class IIa

(Tube style 403 approx. 1200 mm)

For ultrasonic – tube style 403 (incl. silicone sleeve for extra hold) / 1200 mm

2002-02-1-954, Class IIa

(Tube style 403 approx. 300 mm)

For ultrasonic – tube style 403 (incl. silicone sleeve for extra hold) / 300 mm

2002-02-1-932, Class IIa

(Tube style 303 approx. 1200 mm)

For Ultrasonic – tube style 303/1200 mm

2002-02-1-928, Class IIa

(Tube style 303 approx. 300 mm)

For Ultrasonic – tube style 303/300 mm

2002-02-1-914, Class IIa

(Bacteria filter USV2002 disposable)

For Ultrasonic – bacteria filter / disposable item

2002-01-4-000, Class I

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For Ultrasonic – accessory kit for wall variant

2002-01-3-000, Class I

(Accessory kit / USV 2002 table-top model)

For Ultrasonic – accessory kit for table-top variant

2002-00-5-900, Class I

(USV carriage / accessory kit)

For Ultrasonic – accessory kit for carriage variant

2002-00-5-400, Class I

(Carriage USV with divided column)

For Ultrasonic – accessory kit for carriage variant with divided column

2002-00-1-902, Class IIa

Refilling system BioLogic for USV

⚠ ATTENTION!

The use of accessories or components with the USV2010 that are not listed here may result in reduced immunity to interference or increased EMC emissions.

Disposal

11.4 Interface description

All devices or accessories that are combined with the ultrasonic nebuliser must be listed here:

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All devices or accessories must fulfil the conditions of the following interface description under the intended conditions of use according to the intended purpose. Only accessories and consumables suitable for the specific application, which bear the CE mark and which fulfil the characteristics defined in the interface description are permitted. The configuration of the overall system, as well as the verification of its functionality, are the overall responsibility of the authorised personnel.

12. Disposal

At the end of the service life of the product, the ultrasonic nebuliser must be properly disposed of. In doing so, ensure that the materials are carefully separated. The device and accessories do not contain any hazardous goods; the device component materials are fully recyclable.

Whenever the device is passed on, including for disposal, it must be properly disinfected.

This product falls within the scope of EC Directive 2002/96/EC (Waste Electrical and Electronic Equipment Directive). The distributor or a specialist dealer authorised by the distributor will take back the device at the end of its use subject to a charge.





Disposal via household waste collection or the municipal collection points for used electrical equipment is not permitted. For more detailed information regarding legally compliant disposal, please contact your authorised dealer.

Outside the European Economic Area

The applicable national regulations for the disposal and treatment of used equipment must be observed when disposing of this product.

13. Technical data

13.1 Environmental conditions

Storage	
Ambient temperature	 - 20 °C to + 70 °C
Relative humidity, without dewing	 0 % to 90 %
Air pressure	 800 hPa to 1060 hPa
Use	
Ambient temperature	 + 10 °C to + 30 °C
Relative humidity, without dewing	 0 % to 90 %
Air pressure	 800 hPa to 1060 hPa

13.2 Performance data

Voltage	110 V/ 127 V/ 230 V 50/60 Hz
Fuses	T1.0A H 250 V 5x20 mm, 2 units
Power consumption	Ultrasonic nebuliser 50 VA; Tube heating 30 VA
Tube heating resistance	3 Ohm (± 5 %)
Ultrasonic frequency	1.68 MHz ± 5 %
Sound pressure level	approx. 45 dB(A)
Protection class	Protection class I, type B
Protection type	IP 22 (dripping water)
Operation type	Continuous operation, time-controlled
Nebuliser output	max. 3 ml/min
Nebuliser output with sterile water system	up to 2.5 ml/min
Nebuliser chamber volume	0 – “Min” approx. 80 ml; “Min” – “Max” approx. 85 ml; Total volume approx. 170 ml
Droplet size	86% < 5 µm, average 4 µm
Air flow	up to 20l/min with use of the bacteria filter
Device dimensions W x H x D	205 x 190 x 315 mm
Weight	approx. 3.500 kg
Classification according to EC Directive 93/42/EEC Annex IX	Class IIa
UMDNS code	12-719
GMDN	12719

Electromagnetic compatibility

14. Electromagnetic compatibility

14.1 Electromagnetic emissions

Guidelines and manufacturer's declaration - electromagnetic emissions

The device is intended for operation in an environment as specified below. The customer or user of the device should ensure that it is operated in such an environment.

Interference emission	Conformity	Electromagnetic environment – guidelines
HF emissions per CISPR11	Group 1	The device generates HF energy internally and conducts this to the transducer for the purpose of generating mist. The HF emission is very low. However, electronic equipment in the immediate vicinity may be disturbed.
HF emissions per CISPR11	Class B	The device is intended for use in all clinical facilities (qualified personnel).
Emission of harmonics according to IEC 61000-3-2	Class A	
Emitted voltage fluctuations / flicker per IEC 61000-3-3	Is compliant	

14.2 Electromagnetic interference on other equipment

Guidelines for avoiding, detecting and eliminating electromagnetic interference on other equipment

Other electrical / electronic equipment should not be operated in the immediate vicinity of the device or stacked with it. If such an arrangement cannot be avoided, the equipment should be observed as appropriate to ensure it functions correctly. Due to the variety of device functions, the effects can be quite different and sometimes difficult to detect.

Examples are described as follows:

Device

- Type of fault.
 - ☞ Remedy.

Radio, TV appliances

- Noise / crackling in the sound, horizontal stripes in the image.
 - ☞ Increase distance.
 - ☞ Change arrangement.
 - ☞ Change orientation.
 - ☞ Changing reception channel.

Monitoring systems, e.g. Baby listener / cordless telephone

- Noise / crackling in the sound.
 - ☞ Increase distance.
 - ☞ Change arrangement.
 - ☞ Change orientation.
 - ☞ Changing reception channel.

Wireless thermometer / wireless weather stations

- Interference with data transmission, missing or incorrect display.
 - ☞ Increase distance.
 - ☞ Change arrangement.
 - ☞ Change orientation.
 - ☞ Changing reception channel.

General electronic devices

- Malfunction, e.g. stop or change of the intended operating mode.
 - ☞ Increase distance.
 - ☞ Change arrangement.
 - ☞ Change orientation.
 - ☞ Changing reception channel.

14.3 Electromagnetic interference immunity

14.3.1 Guidelines and manufacturer's declarations - electromagnetic interference immunity

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should ensure that it is operated in such an environment.

The environment is described as follows:

Interference immunity

IEC 60601 test level / compliance level

- Electromagnetic environment – guidelines

Discharge of static electricity / ESD per IEC 61000-4-2

- ± 8 kV contact discharge
- ± 15 kV air discharge

- Floors should be made of wood or concrete, or be covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.

Rapid transient electrical disturbances / bursts per IEC 61000-4-5

- ± 2 kV for mains lines
- The quality of the supply voltage should correspond to that of a typical commercial or hospital environment.

Impulse voltages / surges per IEC 61000-4-11

- ± 1 kV series mode voltages
- ± 2 kV common mode voltages
- The quality of the supply voltage should correspond to that of a typical commercial or hospital environment.

Voltage dips, short interruptions and supply voltage fluctuations according to IEC 61000-4-11

- 0 % U_T / 0.5 periods
- 0 % U_T / 1 period
- 70 % U_T / 25 periods
- 0 % U_T / 250 periods

- The quality of the supply voltage should correspond to that of a typical commercial or hospital environment. If the user of the device requires continued operation even in the event of power interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

Remark:

U_T is the mains AC voltage before applying the test level.

Magnetic field with the supply frequency (50/60 Hz) per IEC 61000-4-8

30 A/m

- Magnetic fields at the mains frequency should correspond to the typical values found in the commercial and hospital environment.

Electromagnetic compatibility

14.3.2 Guidelines and manufacturer's declarations - electromagnetic interference immunity, non-life-supporting equipment

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should ensure that it is operated in such an environment.

The environment is described as follows:

Interference immunity

IEC 60601 test level / compliance level

- Electromagnetic environment – guidelines

Conducted HF disturbance variables according to IEC 61000-4-6

6 V_{eff}
150 kHz - 80 MHz

Radiated HF disturbance variables according to IEC 61000-4-3

3 V/m
80 MHz - 2.7 GHz

- Portable and mobile radios should not be used at a distance from the product, including cables, that is less than the recommended safety distance (30 cm) calculated according to the equation applicable to the transmission frequency.
- The field strength of stationary radio transmitters should be lower than the compliance level at all frequencies according to on-site testing.
- Interference may occur in the vicinity of equipment marked with the following symbol:

Remark:

At 80 MHz and 800 MHz the higher frequency range applies.

14.4 Recommended safety distances

Guidelines and manufacturer's declaration / high-frequency wireless communication devices

Test frequency (MHz)	Radio service	Maximum power (W)	Distance (m)	Test level (V/m)
385	TETRA 400 (380...390 MHz)	1.8	0.3	27
450	GMRS 460, FRS 460 (430...470MHz)	5	0.3	28
710, 745, 780	LTE Band 13/ 17 (704...787MHz)	0.2	0.3	9
810, 870, 930	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 (800...960MHz)	2	0.3	28
1720, 1845, 1970	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1/ 3/ 4/ 25, UMTS (1.7...1.99GHz)	2	0.3	28
2450	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band7 (2.4...2.57GHz)	2	0.3	28
5240, 5500, 5785	WLAN 802.11 a/ n (5.1...5.8GHz)	0.2	0.3	9

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All rights reserved and subject to change.
W. Krömker GmbH