Approved
for air travel!

DeVilbissBLU $\approx$ ™

**DV63 DeVilbiss® Standard Plus &
DV64 DeVilbiss AutoPlus® CPAP Series**

CE 0044

EN DeVilbiss BLUE™ Positive Airway Pressure Device Instruction Guide

RX ONLY

Made in USA of US and Imported Parts. Not manufactured with natural rubber latex.

CE DV63E & DV64E series devices are CE0044
DV63A & DV64A series devices are not.

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


Display

PRESSURE RAMP

When activated, this feature gradually ramps up the device to prescription pressure to help you fall asleep easier. Ramp time can be set from 0-45 minutes; refer to Features & Settings for details on adjusting length of time.

To start ramp:

Press RAMP button


Display

Ramp On
(Time Remaining)

To stop ramp: (begin prescription pressure)

Press RAMP button again

Display

Ramp Off

HUMIDIFICATION

Humidification (optional)


Clean, fill and install humidifier chamber daily. **NOTE:** Use distilled water to avoid mineral deposit build up inside the chamber.

To Set Humidification Controls:

1. Select Humidity.

2. Choose humidity setting of OFF or 1 through 5 bars. (NOTE- 5 bars is the highest humidity)

NOTE – For models with heated humidification option, please refer to humidifier guide for device set up.

**if Auto-START/STOP feature is activated

ENGLISH.....	EN - 2
ESPAÑOL	ES - 32
FRANÇAIS	FR - 63

TABLE OF CONTENTS

Quick Start - Operating Your CPAP	EN - 2	Clock Settings Menus	EN - 15
Symbol Definitions	EN - 4	Clock - Current Time	EN - 15
Important Safeguards	EN - 4	Format	EN - 15
Introduction	EN - 6	Wake-Up Alert	EN - 15
Intended Use	EN - 6	Add-Ons	EN - 16
Contraindications	EN - 6	Bluetooth	EN - 16
Set Up	EN - 6	Wireless Modem	EN - 16
Important Parts	EN - 7	Oximeter	EN - 16
Understanding your Display	EN - 7	Other Device	EN - 17
Understanding your Controls	EN - 8	Info	EN - 17
Understanding your Device	EN - 8	Patient Messages	EN - 18
System Assembly	EN - 9	Notifications	EN - 18
Features and Settings	EN - 10	Messages	EN - 18
Main Menu Items	EN - 10	Add-On Device Messages	EN - 18
Quick View Menus	EN - 10	Service Codes	EN - 18
SmartCode® Data Menus	EN - 11	SmartLink® Patient Therapy Management System	EN - 19
SmartCode Rx Submenu	EN - 11	Using an SD Card	EN - 19
CPAP Settings Menus	EN - 12	Travel Information	EN - 20
Rx Pressure	EN - 12	Supplemental Oxygen	EN - 21
Auto Mode	EN - 12	Replacement Items/Spare Parts	EN - 22
CPAP Mode	EN - 12	Cleaning	EN - 23
SmartFlex™ Exhalation Relief Comfort Technology	EN - 12	Maintenance	EN - 25
Ramp Options	EN - 13	Expected Service Life	EN - 25
Humidity Settings	EN - 13	Product Disposal	EN - 25
Tubing Settings	EN - 13	Network/Data Coupling	EN - 25
Master Volume Settings	EN - 14	Troubleshooting	EN - 26
Brightness Settings	EN - 14	Specifications	EN - 28
Enable Menu	EN - 14	Guidance and Manufacturer's Declaration	EN - 30

SYMBOL DEFINITIONS

	It is mandatory to read Operating Instructions prior to use.		Heat
	Consult instructions for use		RTCA / DO-160 Section 21 Category M for battery use only
	Rx Only = Prescription required		Manufacturer
	Class II electrical protection-double insulated		Date of Manufacture
	DC input		Radio Transmitter
	DC Jack		TUV Rheinland C-US approval mark
	Type BF equipment-applied part		European Representative
	Catalog Number		European Rep CE mark
	Serial Number		MR Unsafe – Unsafe for Magnetic Resonance Environment
	IP21		
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU-Waste Electrical & Electronic Equipment		

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING THIS DEVICE.

SAVE THESE INSTRUCTIONS.

The device is to be used only on the instruction of a licensed physician. The system can deliver CPAP therapy or Auto-CPAP therapy. For enhanced pressure relief, the device can also deliver SmartFlex pressure relief. Your home care provider will make the correct pressure settings according to your health care professional's prescription.

When set in the Auto-CPAP therapy, the system will monitor your breathing as you sleep and automatically adjust the pressure to meet your needs. When in CPAP therapy, the system will deliver a continuous, set pressure during the night.

WARNING

A warning indicates the possibility of injury to the user or the operator.

- **Electric Shock Hazard** – Do not use while bathing.
- **Electric Shock Hazard** – Do not immerse this device into water or any other liquid.
- **Electric Shock Hazard** – Do not attempt to open or remove the cabinet; there are no user-serviceable internal components. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.
- Refer to International Standard IEC 60601-1 Ed 3.0 Amendment 1 for safety requirements applicable to Medical Electrical Systems
- Oxygen supports combustion. To avoid possible physical injury, do not smoke while using this device with supplemental oxygen. Do not use this device near hot objects, volatile substances or sources of open flames.
- Always turn on the device before turning on the oxygen source. Turn off the oxygen source before turning off the device. Never allow the oxygen source to run continuously while connected to the device if the device is not in use. If the device is not being used, turn off the oxygen flow.
- CPAP pressure must be set higher than 8 cmH₂O for supplemental oxygen use to prevent patient breathing pattern from forcing oxygen back into the CPAP.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary depending on the pressure setting, patient breathing pattern, mask selection, and leak rate. This warning applies to most types of CPAP devices.
- The device should be used only with vented CPAP masks recommended by DeVilbiss, your physician or respiratory therapist. These include nasal, full face, nasal pillow and nasal cannula masks designed for CPAP use. Do Not use a non-vented mask with this CPAP device.
- To avoid rebreathing of exhaled air, do not use a CPAP mask unless the device is turned on and providing a supply of air. Venting in the mask should never be blocked. When the device is turned on and providing a fresh supply of air, exhaled air is flushed out of the mask vent. However, when the device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most CPAP devices.
- The device is not a life support device and may stop operating with certain device faults or with a power failure. It is intended to be used on spontaneously breathing individuals weighing 66 lbs/30 kg or greater.
- To avoid electric shock, always unplug power cord from wall outlet power source when cleaning or servicing.

- The following surface temperatures may exceed 41 °C under certain conditions:
 - External surface of CPAP 43°C
 - Patient tubing connector 45°C
 - Encoder shaft (with control knob removed) 42°C
 - Humidifier heater plate 65°C
 - External power supply 42°C
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.
- Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.
- The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.
- If you are using a full face mask (a mask covering both your mouth and your nose), the mask must be equipped with a safety (entrainment) valve.
- Do not connect the device to an unregulated or high pressure oxygen source.
- Do not use the device near a source of toxic or harmful vapors.
- Do not use this device if the room temperature is warmer than 40° C (104° F). If the device is used at room temperatures warmer than 40° C (104° F), the temperature of the airflow may exceed 43° C (109° F). This could cause irritation or injury to your airway.
- Do not operate the device in direct sunlight or near a heating appliance because these conditions can increase the temperature of the air coming out of the device.
- Contact your health care professional if symptoms of sleep apnea recur.
- If you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, if it has been dropped or mishandled, if water is spilled into the enclosure, or if the enclosure is broken, disconnect the power cord and discontinue use. Contact your home care provider.
- Periodically inspect electrical cords and cables for damage or signs of wear. Discontinue use and replace if damaged.

CAUTION

A Caution indicates the possibility of damage to the device.

- For proper operation, place system on a flat, sturdy, level surface. Do not place the

device where it can be knocked onto the floor or where the power cord may create a trip hazard.

- Never block air openings of the device. Do not insert objects into any openings or tubes.
- Use only accessories recommended by DeVilbiss.
- The USB port located on the back of the device is for provider use only. The port must only be used with accessories approved for use by DeVilbiss. Do not attempt to attach any other device to this port as it may damage the CPAP or the accessory device.
- Only the DeVilbiss DV6 series Heated Humidifier system is recommended for use with the device. Other humidifier systems may prevent the device from detecting respiratory events and may cause inappropriate pressure levels in the mask.
- The height of the device must be lower than the mask when using a humidifier to prevent water from getting into the mask
- Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in product damage. If this occurs, do not use the device. Disconnect the power cord and contact your equipment provider for service.
- Condensation may damage the device. If this device has been exposed to either very hot or very cold temperatures, allow it to adjust to room temperature (operating temperature) before starting therapy. Do not operate the device outside of the operating temperature range shown in the Specifications.
- Do not place the device directly onto carpet, fabric, or other flammable materials.
- Tobacco smoke may cause tar build-up within the device, which may result in the device malfunctioning.
- Dirty inlet filters may cause high operating temperatures that may affect device performance. Regularly examine the inlet filters as needed for integrity and cleanliness.
- Never install a wet filter into the device. You must ensure sufficient drying time for the cleaned filter.
- Only use a DeVilbiss DC Power Cord and Battery Adapter Cable. Use of any other system may cause damage to the device.
- Oxygen is a prescribed gas and should only be administered under the supervision of a physician.
- The flow setting of the oxygen source must be specified by a physician.
 - Maximum oxygen pressure is 50 psi. Maximum oxygen flow rate is 10 lpm.
- Always turn on CPAP before turning on oxygen flow.
- Always turn off oxygen flow before turning off CPAP.

- Proper filter function is important for the operation of the device and to protect the device from damage.

INTRODUCTION

Intended Use

The DeVilbiss IntelliPAP2/DeVilbiss BLUE Series is intended for use in treating OSA in spontaneously breathing patients 30 Kg (66 lbs) and above by means of application of positive air pressure. The Device is to be used in Home and Healthcare Environments.

Contraindications

Positive airway pressure therapy may be contraindicated in some patients with the following pre-existing conditions:

- severe bullous lung disease
- pneumothorax
- pathologically low blood pressure
- dehydration
- cerebrospinal fluid leak, recent cranial surgery or trauma
- a bypassed upper airway
- hypoventilation

The use of positive airway pressure therapy may be temporarily contraindicated if you exhibit signs of a sinus or middle ear infection. Contact your physician if you have any questions concerning your therapy.

MR Unsafe

Do not bring the device or accessories into in a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the CPAP or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.

Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the CPAP. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.

SET UP

First, unpack your CPAP and identify all items.

Contents

1. Carrying Case
2. Tubing
3. Filter Cover
4. Extra Air Inlet Filters
5. Optional Fine Particle Filter
6. SD Card (optional)
7. Power Supply
8. Power Cord
9. CPAP Device


Instruction Guide (not shown)

NOTE - Contents vary by model. Please refer to the contents list on the device carton for items specific to your model.

NOTE - For models with heated humidification option, please refer to humidifier guide for device set up.


IMPORTANT PARTS

Understanding Your Display



IMPORTANT PARTS

Understanding Your Controls

Understanding Your Device


SYSTEM ASSEMBLY

1. Place system components on a stable surface.

NOTE – For models with heated humidification option, please refer to humidifier guide for device set up.

2. Insert power cord into back of device and wall outlet.

3. Verify air supply port cover is inserted into port on bottom of device.

4. Attach tubing to device.

5. Prepare mask and connect to tubing.

6. Insert SD card into slot (if applicable).

FEATURES AND SETTINGS

Your device has many features and settings to personalize and improve your therapy experience. Many of these are fully adjustable, while others may have limited or no adjustability. The following legend is used throughout this section:

 **Informational Only** - Indicates a non-adjustable setting

 **Enabled/Disabled** – Allows you to turn selected options on or off



 **Locked** - Indicates a setting that can be locked by your provider. Locked  or unlocked  icon will be shown in the upper right side of the Display. Locked settings are only adjustable by your provider.

 **Shown/Hidden** - Indicates a feature that can be hidden by your provider. If a feature is hidden, you will not see it in your device.

 Within menu, select Back to return to previous screen. **NOTE:** Pressing On/Off button while in any menu will exit to therapy screen.




Main Menu Items

Without humidifier option:

Navigating Screens: Simply turn rotary knob to desired item and press knob to select item.

With optional humidifier:


Navigating Screens: Simply turn rotary knob to desired item and press knob to select item.


Quick View Menus

Quick View shows an overview of your therapy, as well as more-detailed information for a preset number of days. There is no user adjustability in Quick View. Your provider will contact you for this information if required by your insurer.

1. Press to select **Quick View**.

2. Overview screen appears.

3. Turn knob to choose number of days to view.


FEATURES AND SETTINGS

☰ SmartCode® Data Menus ⓘ

SmartCode data sets offer a detailed look at your therapy through an encrypted code. The decrypted data includes info required by many insurers to continue coverage of your equipment. Data sets can be decrypted by visiting the www.DeVilbissSmartCode.com website and entering the code. Your provider will contact you for this data set if required by your insurer.

1. Select **SmartCode** option


2. 1-Day data set appears

3. Continue to turn knob to view screens:

- 7 Days
- 30 Days
- 90 Days
- Compliance
- Date
- Serial number
- SmartCode Rx (see following section)


Generating Compliance Reports for Employer/Insurer/Physician

SmartCode reports satisfy most requirements for ongoing compliance/adherence information, and the software retains 2,046 days (5.6 years) of data. We recommend you work with your provider to establish a routine for clearing SmartCode data to ensure current usage information is available. Compliance meter should be cleared only after the required report is generated, submitted, **and approved** by the requestor. Once cleared, a new cycle of 5.6 years begins and the only source for the previous data will be the hard-copy reports.

NOTE Only the provider has access to clear the meter.

SmartCode Rx Submenu

This submenu allows your provider to adjust your prescription through input of an encrypted code. Code may contain letters or numbers and knob may be turned either way to roll through 0-9 and A-Z.

1. Select **SmartCode Rx** to input code supplied by your provider

2. Highlight, then turn knob to select first digit and press.

3. Continue through all digits

4. Write down Verify Code and give to provider


FEATURES AND SETTINGS

CPAP Settings Menus

RX Pressure

Select Rx Pressure

Auto Mode

View your prescription pressure setting

CPAP Mode

View your prescription pressure setting

SmartFlex™ Pressure Relief Comfort Technology

When activated, the SmartFlex pressure relief feature helps make it easier for you to exhale against your prescribed pressure by slightly lowering the pressure during exhalation. If you experience a respiratory event when using the pressure relief function, the device will automatically deactivate this function until normal breathing resumes.

1. Select SmartFlex


2. Select Setting

3. Choose 1-3 or OFF

NOTE – 1 is the least relief and 3 is the most.

4. Press to choose Always On or Ramp Only


Always On = SmartFlex pressure relief is active throughout the night.

Ramp Only = SmartFlex pressure relief is on during ramp only.

5. Select Flow Rounding


6. Choose 0-5 for Inhalation

7. Choose 0-5 for Exhalation

NOTE – 0 is the least rounding and 5 is the most.

(Flow Rounding determines the speed at which the SmartFlex pressure relief changes.)


FEATURES AND SETTINGS

⚙️ CPAP Settings Menus

📶 Ramp Options 🔒 ⓘ

When activated, this feature gradually ramps up the device to prescription pressure to help you fall asleep easier.


1. Select Ramp


4. View Ramp Pressure ⓘ (set by the provider)

2. Choose Ramp Time. (if unlocked).


3. Set Ramp Time 🔒 (From 0-45 minutes in 5 min increments)


🌫️ Humidity Settings

(Setting only applicable if you have the optional heated humidification.)

1. Select Humidity


2. Choose humidity setting of OFF or 1 through 5 bars. **NOTE:** 5 bars is the highest humidity.

🌀 Tubing Settings 🔒



Allows adjustment for optimal therapy based on your tubing dimensions.

1. Select Tubing


2. Select Diameter: 15mm or 22mm

Diameter	15mm	22 mm
Length	1.8 m (6')	1.8 m (6'), 2.4 m (8'), 3.0 m (10'), 3.7 m (12')

3. Select Length if 22 mm Diameter. **NOTE:** These settings must match the tubing you are using for pressure at the mask to be correct.

NOTE: Refer to Replacement Items/Spare Parts section to verify tubing diameter.


NOTE: Bacteria filter option is set by the provider.

CPAP Settings Menus

Master Volume Settings

Allows adjustment of alert volume for patient messages.

1. Select Volume


2. Set Master Volume (sound level) from 0-100% 0 = Off (mute)

Brightness Settings

Allows adjustment of display and controls brightness.


1. Select Brightness

2. Select Display (brightness of display area) or Controls (brightness of top buttons, heater-on indicator, and rotary knob)

3. Select Active or Standby

4. Select Brightness Level

Active = Brightness during user activity (when using buttons or menus). Settings from 10% to 100% (must be equal to or higher than Standby brightness)

Standby = Brightness when no user activity (idle). Settings from 0% to Active (must be equal to or lower than Active brightness)

Display and controls change to Standby brightness after 2 minutes of no activity.


NOTE – If Standby controls brightness is set to 0, the heater on indicator remains on at a low level.

Enable Menu


✓ = Enabled (ACTIVATED)

Allows Auto-Start/Stop feature and Notifications to be activated or deactivated. Activating Auto-Start/Stop allows you to start therapy by breathing into your mask and to stop therapy by removing your mask. Activating Notifications will ensure they are displayed on your screen when generated.


1. Select Enable and then select desired options


2. Auto-STOP **NOTE**- If Auto-STOP is enabled, Auto-START is automatically enabled.

3. Auto-START

4. Notifications (Refer to Patient Messages for details.)


FEATURES AND SETTINGS

Clock Settings Menu


Clock - Current Time

1. Select Current Time


2. Set hour, minutes, and AM/PM (if applicable)

12 24 Format

1. Select 12 or 24 hour Clock Format


Wake-Up Alert


Wake up alert allows your CPAP to also function as an alarm clock.


1. Select Wake-Up Alert


2. Select Wake-Up Alert On

3. Select & Set Alert Time

4. Select Volume

5. Set Wake-Up Alert Volume (sound level) 0=Off (mute)

NOTE – Wake-Up Alert Volume is separate from Master Volume.

If provider hides clock, wake up alert is hidden. If clock is shown, wake up alert can be shown or hidden.

NOTE – If Standby display brightness of 0% (Off) is chosen, the Clock will not be visible when idle. Wake-Up Alert is available and changes the screen to Active upon activation.

FEATURES AND SETTINGS


≡ Add-Ons

This menu allows the connection of DeVilbiss approved devices such as DV6WM Wireless Modem, Nonin® WristOx2® wireless oximeter, and other devices.


✚ Bluetooth®

Bluetooth allows connection of Bluetooth enabled items such as wireless modem or oximeter.


1. Select **Bluetooth** and press to select On.

2. Select search.

3. Searching for Bluetooth devices. Shows list of oximeters in range. Click the item serial number to pair with the CPAP.


For Airline Use – Select **Bluetooth** and press to select Off.


📶 Wireless Modem

Information on optional DV6WM.

1. Select **Wireless Modem**

2. Wireless modem information screen


Shows signal strength of Bluetooth and wireless modem as well as battery status of modem.

NOTE – Refer to DV6WM Instruction Guide for details.


📶 Oximeter

Information on optional Nonin WristOx2 wireless oximeter.

1. Select **Oximeter**


2. Oximeter information screen. **NOTE** – Refer to Oximeter Instruction Guide for details.

NIGHTLY

1. Connect wireless oximeter to CPAP and place probe on finger.
2. Wait until your CPAP display shows the pop-up message "Oximeter Logging Started" which means your CPAP is recording data from the oximeter.

NOTE – If more than one person in the house is using a wireless oximeter, each person must do this check at a different time.


3. After recording is finished, remove finger probe and disconnect the sensor cable from the oximeter. The CPAP will display a message 'Oximeter Logging Stopped' and beep once.

Add-Ons

Other Device

1. Select desired Device.


2. Other Device information screen.




Info

Allows you to see basic device information.

1. Select Info.

2. Turn knob to scroll through Compliance Meter, Serial Number, Firmware, and BIOS screens.

PATIENT MESSAGES

Various messages may be displayed during your therapy; some include an audible alert that is controlled by Master Volume Setting.

1. Notifications

These are friendly reminders to clean or replace components and, if enabled, are generated as follows. Dismiss from Display by Button Press or Knob Movement.

Displayed Notification	Recurrence interval
Air Intake Filter Due for Cleaning	50 hours of use (while breathing)
Mask, Tubing, and Chamber Due for Cleaning	35 hours of use (while breathing)

2. Messages

These identify conditions that may require action by you/your provider or are informational only. Dismiss from Display by Button Press or Knob Movement.

Displayed Message	Condition for Occurrence
Mask Off Please Check Mask Fit	10 seconds after mask off condition detected.
Auto-STOP Detected	When Auto-STOP occurs.
Auto-START Detected	When Auto-START occurs.
Wake-Up Alert <current time>	When Wake-Up Alert activates (alarm clock).
Error with Code, Prescription Not Updated	When SmartCode Rx is entered but incorrect.
Card Detected Please Wait ...	When SD card is detected. Do not remove SD card while this message is displayed.
Card Removed Send to Provider	When SD card is removed after data transfer.
Card Error Contact Provider	If SD card error condition detected.
Settings Not Updated Contact Provider	When prescription change via SD Card fails.
Card Transfer Successful	When valid SD Card inserted and data transfer complete.
Card Detected Ready for use	When valid SD Card is inserted.
Firmware Update Please Wait ...	When SD Card containing firmware update is inserted.

3. Add-On Device Messages

These messages may be displayed when Add-On items are in use, Dismiss from Display by Button Press or Knob Movement.

Displayed Message	Condition for Occurrence
Oximeter Logging Started SpO2 _____ Pulse Rate _____	When Oximeter first detected and logging starts
Oximeter Logging Stopped	When Oximeter signal is lost and logging stops
Oximeter Finger Probe Off	When Oximeter finger probe off condition detected
Oximeter Finger Probe On SpO2 _____ Pulse Rate _____	When Oximeter finger probe off condition corrected
Wireless Modem Low Battery, Please Recharge	When DV6WM is paired and battery level is low

4. Service Codes

These are divided into Critical or Non-Critical. Critical service codes put the unit in a fail-safe (blower off) state.

Displayed Message	Action
Critical Service Code E0X - Where (X) is a number or letter NOTE – Code E01 is always shown in English since language setting is unknown	Contact your equipment provider to return unit for service. Message remains on Display and cannot be dismissed.
Service Code E8X (non-critical - Where (X) is a number or letter	Contact your equipment provider. Dismiss from Display by Button Press or Knob Movement.

NOTE- Service Code alerts are unaffected by notification enable/disable setting.

DEVILBISS SMARTLINK® PATIENT THERAPY MANAGEMENT SYSTEM

DeVilbiss SmartLink technology is included in every device. It is used by your provider in conjunction with our exclusive software to monitor the effectiveness of your therapy and how consistently you use the device. Information can be obtained through the use of an SD card or by direct download to a computer with SmartLink software installed.

Using an SD Card

Your provider may give you an SD card to collect data. It should be installed during device set up and kept in device. Your provider will instruct you on when to return the card.

To install the card

1. Insert card face up. Press until it clicks into place.

2. Message will display and tone will sound.

Card Detected
Ready For Use

Data transfer

1. Device is copying card. Card should remain in device until copying is complete

NOTE – DO NOT remove SD card while showing “Card Detected Please Wait.” SD card may be corrupted if removed during this message.

Card Detected
Please Wait...

To remove the card

1. Press and release card; remove from CPAP.

2. Message will display and tone will sound.

Card Removed
Send to Provider

TRAVEL INFORMATION

NOTE – For details on part numbers for alternate power sources, refer to Replacement Items/Spare Parts.

Air Travel

Your device is RTCA/D0-160 Section 21 Category M approved for use on commercial airlines when using battery power. You must turn off the Bluetooth. Refer to Features & Settings - Add-Ons in the guide for details on Bluetooth.

International Travel

Your device is capable of accepting line voltages of 100-240V~, 50/60Hz. For travel to another country, simply order the correct power cord for that region from your provider or use proper adapter.

Altitude

This device automatically compensates for altitude between 1060 hPa and 700 hPa (~ 1,400 ft below sea level to 9,800 ft). At higher altitudes, breathing patterns change and you may lose some benefit of your CPAP therapy.

Battery/DC Operation

Several DC power source options are available:

1. DeVilbiss offers an optional, custom battery (DV6EB) for use as a backup power source when AC power is not available, when traveling by air, or any time portability is desired. Please refer to the guide included with your battery for details.
2. Your device can be operated from a 12 V DC power source by using the optional DC Cable.
3. Alternately, your device can be powered from a stand-alone battery. Connect to a 12V Deep Cycle Marine Battery using the optional DC power cord and optional DC Battery Clamp-On Adapter.
4. Lastly, instead of connecting 12V DC power directly to the device, you can pass the DC power through an AC Inverter that will deliver AC power to the system. The power rating of the inverter must be at least 200 watts @ 110V~ / 400 watts @ 220V~.

Approximate Battery Operation Run Time

Example Battery Size = 100 W/Hr

	CPAP Only	CPAP with Heated Humidifier		CPAP with Heated Humidifier and PulseDose Module	
		Heater Setting = 3	Heater Setting = 5	Heater Setting = 3	Heater Setting = 5
CPAP Setting (cmH2O)	Minimum Run Time (Hrs)	Minimum Run Time (Hrs)	Minimum Run Time (Hrs)	Minimum Run Time (Hrs)	Minimum Run Time (Hrs)
5	16.9	8.7	5.6	9.5	5.9
10	11.9	7.5	4.7	7.8	5.4
15	8.7	6.4	4.4	7.1	4.7
20	6.8	5.8	4.0	6.1	4.4

NOTE – Run times are at a breath rate of 20 breaths per minute and standard mask leak.

SUPPLEMENTAL OXYGEN

WARNING-OXYGEN USE

- Oxygen supports combustion. To avoid possible physical injury, do not smoke while using this device with supplemental oxygen. Do not use this device near hot objects, volatile substances or sources of open flames.
- Always turn on the device before turning on the oxygen source. Turn off the oxygen source before turning off the device. Never allow the oxygen source to run continuously while connected to the device if the device is not in use. If the device is not being used, turn off the oxygen flow.
- CPAP pressure must be set higher than 8 cmH₂O for supplemental oxygen use to prevent patient breathing pattern from forcing oxygen back into the CPAP.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary depending on the pressure setting, patient breathing pattern, mask selection, and leak rate. This warning applies to most types of CPAP devices.

CAUTION


- Oxygen is a prescribed gas and should only be administered under the supervision of a physician.
- The flow setting of the oxygen source must be specified by a physician.
 - Maximum oxygen pressure is 50 psi. Maximum oxygen flow rate is 10 lpm.

If your physician prescribed supplemental oxygen it can be added one of two ways:

1. Place optional oxygen adapter into air supply port.

2. Insert CPAP tubing into adapter.

3. Attach oxygen tubing to the adapter.

Always turn on CPAP before turning on oxygen flow.

Always turn off oxygen flow before turning off CPAP.

OR

1. Attach oxygen tubing directly to the oxygen port on the mask.

REPLACEMENT ITEMS/SPARE PARTS

Replacement Items/Spare Parts

Air-inlet filter
(4/pk)
#DV51D-602

Fine particle
filter (4/pk)
#DV51D-603

Air supply port
cover
#DV61D-604

Filter Cover
#DV63D-631

Air supply tubing
22 mm x 1.8 m (6') -
#DV51D-629
15 mm x 1.8 m (6') -
#DV61D-629

NOTE: Tubing shown at
right is actual size to easily
verify diameter of your
tubing.

Carrying Case
#DV63D-610

Oxygen Adapter
#7353D-601


Inline Bacteria Filter Package (filter,
elbow, 1.8m (6') 22 mm Tubing) -
#DV51D-631

Optional Heated Humidification

Heated Humidifier with Standard Humidification - #DV6HH

Heated Humidifier with PulseDose® Humidification upgrade - #DV6HHPD

Power Items

Power supply
#DV63D-613

AC power cord
(USA)
#DV51D-606

AC power cord (UK)
#DV51D-608

AC power cord (EU)
#DV51D-607

AC power cord
(Australia)
#DV51D-609

DeVilbiss DV6
Battery
#DV6EB

DC power cord
(cigarette lighter
adapter)
#DV6X-619

DC battery clamp-on
adapter (requires DC
Power Cord DV6X-619)
#DV51D-696

CLEANING

WARNING


To avoid electric shock, always unplug the power cord from the wall outlet power source before cleaning.

CAUTION – Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in product damage. If this occurs, do not use the device. Disconnect the power cord and contact your equipment provider for service.

Device Cabinet, AC Adapter (Power Supply)

Device Cabinet

1. Disconnect from power source.


2. Wipe the cabinet with a clean, damp cloth once per week (7 days).


NOTE – Allow to dry completely before plugging into power source.

AC Adapter (Power Supply)

1. Disconnect from device and power source.

2. Wipe the adapter with a clean, damp cloth every few days.

CLEANING

Filters


Air-Inlet Filter - Check every 10 days and clean as necessary.

CAUTION – Proper filter function is important for the operation of the device and to protect the device from damage.

1. Remove filter cover.

2. Remove the dark outer foam filter from the frame.

3. Wash the filter in a solution of warm water and dishwashing detergent, and rinse with water.

NOTE – If filter is damaged, contact your equipment provider for replacement.

4. Ensure the filter is completely dry before installation and use.

5. Install filter cover on back of device. **NOTE**–If filter cover is missing, device will operate normally.

Optional Fine Particle Filter -

Check every 10 days and replace if dirty. **CAUTION** – Proper filter function is important for the operation of the device and to protect the device from damage.

1. Replace every 30 days.

2. Install the fine particle filter first; then install the standard filter.

NOTE – If the fine particle filter is not installed first according to the directions, the life of the filter will be reduced requiring more frequent replacement.

Tubing

Tubing should be cleaned daily.

1. Remove the tubing from the device and mask.

2. Use a mild detergent (e.g. Dawn® dishwashing liquid) and warm water to clean the inside of the tubing. Rinse the tubing and allow to air dry.

NOTE – If a reduced drying time is desired, connect the tubing to the CPAP device and allow air to circulate through the tubing until the tubing is completely dry. Inspect the tubing to ensure that the tubing is completely dry before use.

NOTE – Replace every 6 months.

Carrying Case

1. Wipe with clean cloth dampened with detergent or disinfectant.

MAINTENANCE

WARNING

Electric shock hazard - Do not attempt to open or remove the cabinet; there are no user-serviceable internal components. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.

Disconnect supply before servicing.

Pressure Accuracy—No routine calibration or service is required provided the device is used in accordance with the manufacturer's directions. **NOTE**—*Some countries require periodic servicing and calibration of this kind of medical device. Contact your equipment provider for further information.*

Standard Air-Inlet Filter—Check the filter every 10 days and clean it if necessary. Replace the filter every 6 months or sooner if it is damaged. The standard filter is designed to filter particles larger than 3.0 microns in size. The device **MUST** have the standard filter in place during operation.

Optional Fine Particle Filter—Check the filter every 10 days and replace it if it is dirty or damaged; otherwise replace the filter every 30 days. The optional fine particle filter is designed to filter particles as small as 0.3 microns in size.

Tubing—Check the tubing daily. Replace every 6 months.

EXPECTED SERVICE LIFE

- CPAP - 5 years
- CPAP Tubing - 6 months

PRODUCT DISPOSAL

The outer packaging is made of environmentally friendly materials that can be used as secondary raw materials. If you no longer need this packaging, take it to your local recycling and waste disposal facility according to the applicable regulations for your location.

The equipment, including accessories and internal components, does not belong in your regular household waste; such equipment is manufactured from high-grade materials and can be recycled and reused. The internal motor, speaker, PC boards, display board, and wiring should be removed from the unit and recycled as electronic waste. The remaining plastic components should be recycled as plastics.

The European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) requires that electrical and electronic equipment be collected and disposed of separately from other unsorted municipal waste with the aim of recycling it. The crossed out waste bin symbol indicates that separate collection is required.

NETWORK/DATA COUPLING

- Connection of the DV6 series CPAP and accessories to a network/data coupling that includes other equipment could result in previously unidentified risks to patients and operators, and the responsible organization should identify, analyze, and control such risks;
- Subsequent changes to network/data coupling may introduce new risks, requiring new analysis;

Changes to network/data coupling include:

- Changes in network/data coupling configuration
- Connection of additional items to network/data coupling
- Disconnecting items from network/data coupling
- Update of equipment connected to network/data coupling
- Upgrade of equipment connected to network/data coupling

TROUBLESHOOTING

WARNING

Electric Shock Hazard - Do not attempt to open or remove the cabinet; there are no user-serviceable internal components. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.

Issue	Possible Cause	Remedy
Nothing shows on the display.	1. Device is not plugged in or the power cord is not fully inserted.	1a. Verify that power cord is firmly connected to the device and the power outlet. 1b. If you are using a DC power source, make sure the cable connections are secure. Verify that battery is charged.
	2. Optional Humidifier-The flow generator is not fully seated on humidifier cradle.	2. Refer to assembly instructions to ensure proper contact is made.
	3. There is no outlet power.	3. Locate a suitable functioning power source.
The device does not start when breathing into the mask.	1. The Auto-START feature has been disabled.	1a. Use the ON/OFF button to start and stop device. 1b. If shown, go to CPAP setting then Enable menu.
	2. There is no power to the device.	2. Verify that the power cord is firmly connected to the device and the power outlet.
	3. Breath is not deep enough for the Auto-ON feature to detect.	3. Breathe deeply in and out to start the device.
	4. You are using a full face mask that has an anti-asphyxia valve.	4. Auto-START may not work because your exhaled breath escapes through the open valve. Use the ON/OFF button to start and stop the device.
	5. The air supply port plug is missing or not fully inserted.	5. Ensure the air supply port plug is fully inserted into the device.
	6. The optional humidifier chamber is not fully engaged or is missing.	6. Slide the humidifier chamber fully into cradle. Verify the release latch clicks into place.
The airflow has stopped unexpectedly during use or reports mask leak.	1. The Auto-STOP feature has detected large airflow because of loosely fitting mask.	1. Ensure you have a good mask fit; adjust mask and headgear.
	2. During use, your mouth opens and you begin to mouth-breathe.	2. Contact equipment provider for a chin strap or different mask to prevent mouth-breathing.
	3. The optional humidifier chamber is not fully engaged or is missing.	3. Slide the humidifier chamber fully into cradle. Verify the release latch clicks into place.
	4. The air supply port plug is missing or not fully inserted.	4. Ensure the air supply port plug is fully inserted into the device.
The display shows a service code.	An error has occurred in the device and requires servicing.	Contact your equipment provider for service.
Apnea symptoms have recurred.	1. Air filter may be dirty.	1. Clean or replace the air filter and relocate device away from drapes or other dusty surfaces.
	2. Apnea condition has changed.	2. Contact your physician or equipment provider.

Issue	Possible Cause	Remedy
Skin becomes irritated where mask contacts face.	1. Headgear is too tight or improperly adjusted.	1. Loosen headgear to reduce contact pressure on face.
	2. Your mask may not be sized properly or is not the most appropriate shape for you.	2. Contact physician or equipment provider.
Dryness of throat or nose.	Inadequate humidity.	1a. Add a DeVilbiss DV6 heated humidifier. 1b. Increase the heater setting on the optional humidifier.
Water condensation collecting in the hose causing a gurgling sound when using optional humidifier.	1. Humidification is too high.	1. Reduce humidifier heater setting.
	2. Room temperature fluctuates from higher to lower levels through the night.	2. Increase room temperature.
Air from flow generator seems to be too warm.	1. Air filters are dirty.	1. Clean filters.
	2. Air inlet port is blocked.	2. Unblock air inlet.
	3. Room temperature is too high.	3. Lower room temperature.
	4. Device is located near a heat source.	4. Locate device away from heat source.
	5. Optional humidifier heater setting is too high.	5. Lower heater setting.
Nasal, sinus or ear pain, runny nose	You may have a reaction to air flow pressure.	Discontinue use and contact physician.
Heater setting is not visible on display	1. The flow generator is not fully seated onto the humidifier cradle.	1. Refer to assembly instructions to ensure proper contact is made.
	2. Your CPAP device is not equipped with the optional humidifier.	2. Add a DeVilbiss DV6 series heated humidifier

SPECIFICATIONS

CPAP

Size.....	3.7" (9.4 cm) H x 6.1" (15.5 cm) W x 5.9" (15 cm) D
Weight.....	1.93 lbs (0.9 kg) CPAP only; 3.8 lbs. (1.75 kg) including humidifier
Electrical Requirements AC.....	100-240V~, 50/60 Hz
Electrical Requirements DC.....	11-17 VDC, 5.2 Amps
Maximum Power Consumption (including humidifier).....	65 watts max from AC power source
Typical Power Consumption with Humidifier.....	25 watts
Typical Power Consumption without Humidifier.....	10 watts
Pressure Range.....	3-20 cmH ₂ O
Operating Temperature Range.....	41°F to 104°F (5°C to 40°C)
Operating Humidity Range.....	15% to 93% RH non-condensing
Operating Atmospheric Conditions.....	1060 hPa - 700 hPa (~ 1,400 ft below sea level to 9,800 ft)
Storage & Transportation Temperature Range.....	-13°F to +158°F (-25°C to 70°C)
Storage & Transportation Humidity Range.....	15% to 93% RH non-condensing
Maximum Limited Pressure.....	30 cmH ₂ O under normal use
Sound Pressure Level (tested per ISO 17510-1:2007).....	26.6 dBA
Sound Power Level.....	34.6 dBA

Wireless Specifications:

This medical device contains a radio transmitter. The Bluetooth radio integrated into this device is active when the Bluetooth icon () on the LCD display is illuminated. Bluetooth is used to connect your medical device wirelessly to approved accessories. Turn off the Bluetooth radio when not using a wireless accessory, see ADD ONS section of this guide. If you notice any unexplained changes in the performance of the wireless function or your device, turn off the Bluetooth radio to see if that is the cause. You may experience issues with wireless technology such as delays in displayed or stored data, artifact or gaps in the data due to information lost during wireless transfer.

Radio Technology.....	Bluetooth 2.1 +EDR and Bluetooth 4.0
Bluetooth Power Class.....	1.5
Network Topology.....	Point to Point
Bluetooth Profile Supported.....	SPP
Effective Range.....	50 meters (line of sight)
Effective Radiated Power.....	10 dBm (100 mW)
Radio frequency band (Tx and Rx).....	2.402 GHz ~ 2.480 GHz
Minimum Separation Distance (to other RF transmitters).....	1 cm (0.4 in.)

Security Requirements:

Authentication.....	Enforced on all data channels (outgoing and incoming)
---------------------	---

Filter Specifications

Standard Filter > 3.0 micron particles
 Optional Fine Particle Filter > 0.3 micron particles

Dynamic Short-Term Pressure Accuracy with and without DV6HH Humidifier (Per ISO 17510-1:2007)

Waveform	Volume [mL]	Breath Rate [min-1]	Set Pressure [cmH ₂ O]	Pressure Accuracy CPAP only [Pk-Pk cmH ₂ O]	Pressure Accuracy CPAP with Humidifier [Pk-Pk cmH ₂ O]
Sinusoidal cycle with I:E = 1:1	500	10, 15, 20	3.0	0.5	0.5
			7.0	0.5	0.5
			12.0	1.0	1.0
			16.0	1.0	1.0
			20.0	1.0	1.0

Maximum Flow-Rate (per ISO 17510-1:2007)

Test Pressures					
	3.0 cmH ₂ O	7.0 cmH ₂ O	12.0 cmH ₂ O	16.0 cmH ₂ O	20.0 cmH ₂ O
Measured pressure at the patient connection port (cmH ₂ O)	2.9	6.6	11.4	15.0	19.0
Average flow at the patient connection port (l/min)	78.2	123.5	166.8	157.7	138.0

Long Term Static Pressure Accuracy +/- 0.5 cmH₂O

Additional Specifications

Equipment classification with respect to protection from electric shock Class II

Degree of protection from electric shock Type BF Applied Part

Degree of protection against ingress of liquids IP21 Ingress Protection – Protected against finger access to hazardous parts;
 protected against vertically falling water drops

Mode of operation Continuous

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. See instructions for use with supplemental oxygen.

US Patents

SmartCode Technology is protected under US Patent 8649510

Approvals

Meets RTCA/DO-160 DO-160 Section 21 Category M for battery operation only Airline Use

DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems

This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Enforcement – Guidance	
RF Emissions CISPR 11	Group 2	The DeVilbiss DV63 and DV64 Series CPAP must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.	
RF Emissions CISPR 11	Class B Radiated and Conducted Emissions	The DeVilbiss DV63 and DV64 Series CPAP is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonics IEC 61000-3-2	Class A		
Flicker IEC 61000-3-3	Complies		
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8kV contact ±15kV air	±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2kV on AC Mains	±2kV on AC Mains	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz	V1 = 3 Vrms 6 Vrms on ISM & Amateur Bands	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: $D = (0.4)\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	E1 = 10 V/m	$D = (0.4)\sqrt{P}$ 80 to 800 MHz $D = (0.7)\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment marked with the following symbol: 

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting

This device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz outside ISM bands $D = (0.4)\sqrt{P}$	80 to 800MHz $D = (0.4)\sqrt{P}$	800 MHz to 2.5 GHz $D = (0.7)\sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.1	1.1	2.2
100	3.5	3.5	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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13/17, Rue Joseph Priestley • 37100 Tours
FRANCE
+33 (0) 2 47 42 99 42

DeVilbiss Healthcare Pty. Limited

15 Carrington Road, Unit 8
Castle Hill NSW 2154
AUSTRALIA
+61-2-9899-3144

DeVilbiss Healthcare GmbH

Kamenzer Straße 3 • 68309 Mannheim
GERMANY
+49 (0) 621-178-98-0

SE-DV64-1 Rev B

DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA
800-338-1988 • 814-443-4881 • www.DeVilbissHealthcare.com

DeVilbiss®, IntelliPAP®, AutoAdjust®, AutoPlus®, PulseDose®, SmartLink®, and SmartCode®, are registered trademarks of DeVilbiss Healthcare. DeVilbiss BLUE™, SmartFlex™ and SmartCode Adherence Score™ are trademarks of DeVilbiss Healthcare. Bluetooth® is a registered trademark of Bluetooth SIG, Inc. Nonin® and WristOx2® are registered trademarks of Nonin Medical.

© 2015 DeVilbiss Healthcare LLC. 10.15 All Rights Reserved.

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DeVilbiss Healthcare Ltd
DeVilbiss Healthcare SAS