
Bambus

ROLLSTUHLWAAGE / WHEELCHAIR SCALE

Art.Nr. 900 400 100

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




Instructions for Use


Please read these instructions carefully before operating the device.

Use the correct connector and the correct voltage: 110V or 220V!!!

Do not clean the display unit under running water.


Do not open the display unit. In case of problems, please contact Customer Service or your dealer.


Do not touch the keypad with sharp objects.

Startup

1. Remove transport locks between the measuring cells and the wooden plate

2. Attaching the ramp (see photo)

3. Put the 5 rubber feet on the ground, they must all be evenly on the ground

The scale is powered either by an integrated battery or by a power charger (input 230V ~50 Hz, output 12V = 100mA).

- Place the scale on a firm, even surface.
- Connect the platform of the scale with the socket plug provided for this purpose and with the display unit.
- For mains operation - Connect the display unit via the connection jack provided for this purpose with the power charger connected to the power supply.
- Turn on the scale using the on/off switch on the back of the display unit.
- The scale is now ready for use.

Warning: Only use the original power charger provided. The power charger is not waterproof and must therefore not come into contact with water. If the power charger becomes wet or has visible damage, the scale should no longer be operated with the power charger. In this case, contact your dealer or service partner.

Your wheelchair scale is delivered ready to use. You can use it immediately. If the weight is not displayed correctly, you can repeat the configuration and/or calibration.

Weighing with Wheelchair

The memory function allows you to save a fixed value which will be subtracted during a weighing procedure.

- Use the on/off switch to turn the scale on. Following a function test, [0.0] appears in the display in kg mode.
- Roll the empty wheelchair onto the scale.
- Save the wheelchair weight determined by pressing the →0← button; the scale is also set to zero.
- Roll the wheelchair off the scale.
- The user can now roll the wheelchair onto the platform.
- The user's weight is shown on the display and the saved weight was automatically subtracted.

Select Display Mode

To change the display mode, hold down the (*) button for 2 seconds. The different display modes are as follows:

In normal weighing mode:

Hold down the (*) button for 2 seconds: the weight unit changes to pounds (lb).

Hold down the (*) button once again for 2 seconds; the display is temporarily set to a more detailed display mode "X10". (A more exact weight is determined (for example: two decimal places)).

Hold down the (*) button for 2 seconds, display number of cumulations [n 3]

Hold down the (*) button for 2 seconds, display the upper 4 digits of the cumulative weight [H0]

Hold down the (*) button for 2 seconds, display the lower 4 digits of the cumulative weight [L 1085]

Hold down the (*) button for 2 seconds, [AUt 0], cumulation mode is displayed.

Hold down the (*) button for 2 seconds, [bt 95] battery capacity is displayed.

Hold down the (*) button for 2 seconds, return to regular weighing mode.

Zero Display

When the weight is stable, hold down the 0 button for 2 seconds; [0.00] is displayed and the zero state LED lights up.

Tare Weight and Deleting Tare Weight

Tare weight: When the tare status LED is off and the weight is stable, press the →0← button to save the tare weight and switch to net weight mode; the tare

status LED lights up.

Delete tare weight: When the tare status LED lights up, press the **→0←** button; the mode switches to gross weight mode and the tare weight is deleted. The tare status LED is off.

Energy-Saving Mode (AUtP=X1)

If the device is not used for more than 30 minutes, the display brightness automatically decreases and the display [_] appears.

When a button is pushed or the weight on the platform changes, the display unit automatically returns to regular weighing mode.

Automatic Storage

Press and hold the (*) button six times, [AUt-0] set cumulation mode.

[AUt - 0]: Automatic cumulation is switched off.

[AUt - 1]: Automatic cumulation is switched on and cumulates when weight is added.

[AUt - 2]: Automatic cumulation is switched on when the weight is removed, if you wish to set the automatic cumulation mode.

Delete Storage

Press and hold the (*) button three times; [n 0] is displayed.

Press the **→0←** button, delete cumulation data.

Battery Warning Signal

If the battery has less than 20% of the full charge left, the battery warning display [_] appears. Then [bAtt] appears and the display turns itself off after approximately 2 hours. Connect the power charger to charge the battery.

Switching the Weight Display

Either kilograms (kg) or pounds (lb) can be selected as needed for the weight display.

The scale is factory-set to kg mode.

To do this, hold down the (*) button for 2 seconds. Weight unit is pounds (lb), in order to change the weight unit to pounds.

After every changeover, check whether you are in the desired weighing mode.

Restore Factory Settings

While the device is switched off, press the **0** button and switch on. The display performs a segment test. Now perform the following steps:

Press the **→0←** button; [CAL SP] is displayed.

Press the **→0←** button; [-SEt-] is displayed.

Press the **→0←** button; [-A-d] is displayed.

Press the **→0←** button; [-FACT-] Restore factory settings.

Factory settings: d=0.01 / FS = 150.00 / FLt=0 / AUtP=10 / ADr=00 / b=2400 / AUt=0

Delete storage: (n=0 / H=0 / L=0)

Press the (*) button to confirm [0.0] is displayed. You are back in normal weighing mode.

Maintenance and Cleaning

Before cleaning, disconnect the device from the supply voltage. Do not use any abrasive cleaning agents for cleaning. Only use a moistened cloth or conventional disinfectants. Ensure that no liquids penetrate the device!

Attention: Your scale is equipped with a 6V battery and a 230 / = 12V power charger.

- Always connect the scale to the power charger BEFORE you connect the power charger to the power supply.
- Disconnect the power charger from the power supply BEFORE you disconnect the power supply unit from the scale.

Safety Instructions

Safe operation is no longer guaranteed:

- If the power supply unit has visible damage
- If the power supply unit is no longer working
- After being stored for longer periods of time in humid areas

For your own safety, inform your dealer or service partner in these cases.

Error Message

[_XXXX] Appears when the battery charge level is less than 20%.

[bAtt] Appears as a battery charge level warning: the battery must be charged.

[OUEr] appears when the weight is greater than the capacity final value + 9d.

[-OUEr] appears when the weight is less than -2% of the capacity final value.

[Fehler] appears when the calibration is incorrect. The scale should be recalibrated.

Storage and Transport Conditions

Keep all components of the packaging in the event a return is necessary in order to avoid any potential damage during transport.

Disconnect all connected cables before shipment. This helps to avoid unnecessary damage.

Guarantee

This product will be repaired or replaced if proof is provided that failure occurred within 2 years from the date of purchase (sales receipt is required) due to poor workmanship or defective material. This does not apply to any moving parts such as, for example, batteries, cable, power supply units, rechargeable batteries, etc. This guarantee does not include any abnormal wear and tear or damage caused by an accident or misuse. The guarantee is not valid for devices opened by unauthorized persons. The customer's rights as regulated by law are not affected by this agreement. The dealer from whom the scale was purchased is responsible for claims, customer service and replacement parts.

CE Marking

CE Compliance

This device is interference-free in accordance with the valid EC directive 2004/108/EC.

Note: Under extreme electromagnetic influences, e.g. when operating radio equipment in the direct vicinity of the scale, the displayed value may be affected. After the interference has ended, the product can once again be used as intended; it may need to be switched back on, where required.

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Battery - Disposal

Batteries should not be disposed of with household waste. As a consumer, you are legally obligated to return used batteries. You may turn in your old batteries at public collection points in your community or anywhere where batteries of the relevant type are sold.

Disposal of used electric and electronic devices:

The symbol on the product or its packaging indicates that this product cannot be treated like regular household waste but instead it must be brought to a collection point for recycling electric and electronic devices. Further information can be obtained from your community, the municipal waste disposal service or the store where you purchased the product.

Technical Data

Max. capacity: 300 kg / 660 lbs

Division: 100 g / 0.2 lb

Weighing range: 1 kg – 300 kg / 2.0 lbs – 660 lbs

Display: Easy-to-read 25 mm LCD display

Dimensions: 800 x 800 mm

Weighing platform: Wood with nonskid plastic coating

Power supply: 6V battery and power charger

Zero point tracking: Automatic


Operating temperature: +10° to +35°C

Weight: approx. 11.5 kg

A/D resolution: 1,000,000, A/D sampling rate: 40 times/second, actuating voltage: DC 5V; up to 4x 350 Ω load cells or 8x

700 Ω load cells, high-definition (Å~10) display mode, 6-bit display: 0.8-inch LED selectable display resolution.

Type Label

1. Product name
2. Article number
3. Serial number