

Bath seat Bobby

Item No. 12488KD

drive

DeVilbiss
HEALTHCARE



- The ideal bathing aid for narrow bathtubs
- For placing on the edge of the bathtub
- Makes bathing safer and more comfortable
- Powder-coated aluminum frame
- Fits in all standard bathtubs
- Comfortable, non-slip seat surface with hygienic cutaway
- Mounting arms with non-slip rubber caps
- Available exclusively with backrest
- Easy to clean

Technical Data

Order number	12488KD
Colour	white
Overall width	67 - 78 cm
Overall depth	46 cm
Seat width	41 cm
Seat depth	41 cm

Adjustment (width of bathtub)	45.5 - 57.5 cm
Backrest height	36 cm
Total weight	3 kg
Max. Load	120 kg
PU Pal.	32 Pal.
PU pc.	4 Stk.