

Swivel aid Vitaturn K

Item No. 412100100



- Positioning and rotating made easy
- Light and mobile seat surface provides comfort
- Safety thanks to non-slip surface
- The ultimate bath aid – thanks to the unparalleled height of just 15 mm
- Vitaturn meets the highest quality standards
- The skin-friendly coating is easy to clean
- Can be cleaned in the dishwasher at up to 60 °C

Technical Data

Order number	412100100
Colour	gray
Overall height	1.5 cm
Seat plate diameter	32,2 cm
Base width	32,3 cm

Base length	32,3 cm
Total weight	700 g
Overall depth	0 cm
Max. Load	140 kg