

Swivel and transfer aid Vitaturm L

Item No. 412100200

drive

DeVilbiss
HEALTHCARE



- Simplifies transfer, positioning, and rotating
- Safety thanks to non-slip surface
- The ultimate bath aid – thanks to the unparalleled height of just 19 mm
- Vitaturm meets the highest quality standards
- The skin-friendly coating is easy to clean
- Can be cleaned in the dishwasher at up to 60 °C

Technical Data

Order number	412100200
Colour	white-gray
Overall height	1.9 cm
Seat plate diameter	32,2 cm
Base width	30 cm
Base length	52 cm

Total weight	1.2 kg
Overall depth	0 cm
Max. Load	140 kg
PU Pal.	200 Pal.
PU pc.	10 Stk.