

Raised toilet seat TSE 150

Item No. 5201000



- Two hygienic cutaways and integrated seat
- Recessed grip on the lid makes it easy to open
- Folding, non-slip armrests
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces simplify cleaning
- Delivery includes assembly gloves
- Delivery in practical take-away packaging

Technical Data

Order number	5201000
Colour	white
Height adjustment	6.5 / 10 / 13.5 cm
Overall width	56 cm
Overall length	48.5 cm
Overall height	41 cm
Opening width	23 cm

Opening length	25 cm
Width between armrest	48 cm
Total weight	5.7 kg
Max. Load	150 kg
PU Pal.	32 Pal.
PU pc.	4 Stk.