

Commode TS 130

Item No. 530100100



- Comfortable, safe, and particularly stable
- Large, soft seat surface and backrest for increased comfort
- Seat and cover can be folded up and removed
- Standard receptacle in commercially available design
- Six-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Upholstery seams sealed underneath
- Padded armrests
- Easy to clean

Technical Data

Order number	530100100
Overall width	58 / 62 cm
Overall length	57 cm
Overall height	83,5 - 98,5 cm
Seat width	43 cm
Seat depth	46 cm
Seat height	50 - 64 cm

Seat height adjustment	6
Frame height	35 - 50 cm
Backrest height	36 cm
Total weight	8.5 kg
Max. Load	130 kg
PU Pal.	16 Pal.
PU pc.	2 Stk.

Commode TS 130

Item No. 530100100

drive

DeVilbiss
HEALTHCARE

Accessories



Bucket with lid

Item No.530100090



Seat for TS 130

Item No.530100112



Seat cover for TSE 130

Item No.530100120