

Shower wheelchair/wheeled commode DuoMotion 24 XL

Item No. 530500400



- All the advantages of the standard model
- Additional cross strut ensures maximum stability
- With large 24" drive wheels and hand rims
- Increased mobility and independence
- Including 4 x 5" wheels
- Standard soft seat opens at the front
- Including cover for the hygienic cutaway
- Top quality thanks to powder-coated aluminum
- Standard receptacle in commercially available design
- Machine-washable, adaptable backrest with quick fasteners
- Four-position adjustable seat height

Technical Data

Order number	530500400
Colour	gray
Overall width	73 cm
Overall length	109 cm
Overall height	95 - 102.5 cm
Seat width	56 cm
Seat depth	45 cm
Seat height	51 - 58.5 cm

Seat height adjustment	4
Backrest height	36.5 cm
Clearance height	41.5 - 49 cm
Footrest adjustment	32 - 42 cm
Total weight	22.8 kg
Max. Load	200 kg
PU Pal.	6 Pal.

Shower wheelchair/wheeled commode DuoMotion 24 XL

Item No. 530500400

drive

DeVilbiss
HEALTHCARE

Accessories



Bucket with lid

Item No.530100090



Headrest for DuoMotion & DuoMotion XL

Item No.530500318