

Shower wheelchair/wheeled commode DuoMotion XL

Item No. 530500500



- Additional cross strut ensures maximum stability
- Standard soft seat opens at the front
- Including cover for the hygienic cutaway
- Top quality thanks to powder-coated aluminum
- Standard receptacle in commercially available design
- Machine-washable, adaptable backrest with quick fasteners
- Folding armrests
- The footrests with heel straps can be folded, swiveled, and removed
- Four-position adjustable seat height

Technical Data

Order number	530500500
Colour	gray
Overall width	67.5 cm
Overall length	102 cm
Overall height	95 - 102.5 cm
Seat width	56 cm
Seat depth	45 cm
Seat height	50 - 57.5 cm

Seat height adjustment	4
Backrest height	36.5 cm
Clearance height	41.5 - 49 cm
Footrest adjustment	39 - 49 cm
Total weight	16.8 kg
Max. Load	200 kg
PU Pal.	6 Pal.


Shower wheelchair/wheeled commode DuoMotion XL

Item No. 530500500

drive

DeVilbiss
HEALTHCARE

Accessories



24 drive wheel retrofit set

Item No.530500300



Bucket with lid

Item No.530100090



Headrest for DuoMotion & DuoMotion XL

Item No.530500318