

Shower stool Delphi Turn

Item No. 540400000



- Seat surface swivels easily
- Makes transfer and daily hygiene easier
- Skin-friendly coating
- Three-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Diameter of swiveling seat surface: 32.5 cm
- Anodized aluminum tube legs
- Easy to clean

Technical Data

Order number	540400000
Colour	white-gray
Overall width	48 cm
Overall length	48 cm
Overall height	56.2 cm
Seat plate diameter	32,5 cm

Seat height	46.2 - 56.2 cm
Seat height adjustment	3
Total weight	3.1 kg
Max. Load	140 kg
PU Pal.	80 Pal.
PU pc.	10 Stk.