

# Bath bench Benny

Item No. 550100100

**drive**

**DeVilbiss**  
HEALTHCARE



- Slip-resistant seat surface with drainage holes
- Makes getting into the bathtub easy and safe
- Easily adjustable to fit the bathtub width
- Incl. grab handle and integrated soap dish
- Simple and fast installation without tools
- Rubber feet ensure safety and protect against scratches
- Easy to clean
- Delivery in practical, customer-friendly take-away packaging

## Technical Data

|                               |            |
|-------------------------------|------------|
| Order number                  | 550100100  |
| Colour                        | white      |
| Overall width                 | 35 cm      |
| Overall length                | 68 cm      |
| Adjustment (width of bathtub) | 40 - 63 cm |

|              |         |
|--------------|---------|
| Total weight | 2.1 kg  |
| Max. Load    | 150 kg  |
| PU Pal.      | 75 Pal. |
| PU pc.       | 5 Stk.  |