

Shower chair DSR 130

Item No. 550200300



- Comfortable, safe, and stable
- Large seat surface with drainage holes
- Hygienic cutaway simplifies personal hygiene
- With armrest and backrest
- Five-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Powder-coated aluminum frame
- Easy to clean

Technical Data

Order number	550200300
Colour	white
Overall width	54 - 55 cm
Overall height	86 - 96 cm
Overall depth	55 - 57 cm
Seat width	41 cm
Seat depth	41 cm

Seat height	42 - 52 cm
Seat height adjustment	5
Backrest height	48 cm
Total weight	3.48 kg
Max. Load	130 kg
PU Pal.	16 Pal.
PU pc.	2 Stk.