

Shower stool Delta 150

Item No. 550300100



- Very light and space-saving
- Safe and stable
- Simple, tool-free assembly
- Structured surface ensures secure hold
- Seven-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum profile legs
- Drainage openings
- Easy to clean

Technical Data

Order number	550300100
Colour	white
Overall width	44 cm
Overall length	44 cm
Overall height	59 cm
Seat plate diameter	37 cm

Seat height	44 - 59 cm
Seat height adjustment	7
Total weight	1.7 kg
Max. Load	150 kg
PU Pal.	96 Pal.
PU pc.	4 Stk.