

Swiveling bath seat BES 120

Item No. 550600200



- Makes getting in and out of the bath easier and safer
- Easily rotatable 360° in four steps
- Automatic locking at 90°
- Fits in all standard bathtubs
- Fixation using setscrews with end caps
- Comfortable, non-slip seat surface with hygienic cutaway
- Mounting arms with non-slip plastic sheathing
- Powder-coated aluminum frame
- Easy to clean

Technical Data

Order number	550600200
Colour	white
Overall width	74 - 79 cm
Overall height	59 cm
Overall depth	40 cm
Seat width	40 cm
Seat depth	40 cm

Seat adjustment	4 x 90°
Adjustment (width of bathtub)	58 - 63 cm
Backrest height	46 cm
Total weight	4.3 kg
Max. Load	120 kg
PU Pal.	24 Pal.
PU pc.	2 Stk.