

Bath support base frame

Item No. 555200100



- Height adjustment for a comfortable position in the bathtub that is easy on the back
- Suitable for all three sizes
- Adjustable to two height positions (12.5 cm and 27 cm)
- Easy-care plastic
- Low weight
- Space-saving stowage thanks to folding legs
- Click system enables easy set-up and dismantling

Technical Data

Order number	555200100
Colour	white-green
Overall width	44 cm
Overall length	78 - 98,5 cm

Overall height	15,8 - 28,8 cm
Weight	4.1 kg
Max. Load	85 kg
PU Pal.	20 Pal.