

Walking frame Gecko rigid

Item No. 700100100

drive

DeVilbiss
HEALTHCARE



- Rigid, stable aluminum walking frame
- Ideal for safe walking
- Space-saving folding at the touch of a button
- Height adjustment enables optimal adaptation to the user
- No slipping thanks to rubber caps on the legs
- Low weight
- Can be used over the toilet bowl

Technical Data

Order number	700100100
Colour	silver
Overall width	58 cm
Overall length	47 cm
Overall height	81 - 98.5 cm
Width between handles	42,5 cm
Width between rear legs	52 cm
Width between front legs	48 cm

Folded width	52 cm
Folded height	84 cm
Folded depth	12 cm
Total weight	2.4 kg
Max. Load	120 kg
PU Pal.	48 Pal.
PU pc.	4 Stk.