

Walking frame Gecko reciprocal

Item No. 700100200

drive

DeVilbiss
HEALTHCARE



- Reciprocal aluminum walking frame
- Thanks to joints, the walking frame is flexible and offers the user two supporting legs
- Ideal for safe walking training
- Space-saving folding at the touch of a button
- Height adjustment enables optimal adaptation to the user
- No slipping thanks to rubber caps on the legs
- Low weight
- Can be used over the toilet bowl

Technical Data

Order number	700100200
Colour	silver
Overall width	52 cm
Overall length	45.5 cm
Overall height	79,5 - 97 cm
Width between handles	43 cm
Width between rear legs	52 cm
Width between front legs	49 cm

Folded width	54 cm
Folded height	83 cm
Folded depth	12 cm
Total weight	2.7 kg
Max. Load	100 kg
PU Pal.	48 Pal.
PU pc.	4 Stk.