

Transport chair Expedition Plus

Item No. EXP002BL

drive

DeVilbiss
HEALTHCARE



- Light, versatile, and folding
- Low weight makes transport easier
- Including tilt aid
- Dual-function brake (for braking and parking)
- Backrest can be folded away for space-saving stowage
- Cushioned backrest and seat surface made of washable nylon
- Leg supports can be swiveled aside and removed
- Foot plates can be folded up
- Extra-large 8" steering wheels
- Including lap belt

Technical Data

Order number	EXP002BL
Colour	black-blue
Overall width	56 cm
Overall length	93 cm
Overall height	96 cm
Seat width	42 cm
Seat depth	42 cm
Seat height	- 49 cm

Steering wheels	8"
Folded height	75,5 cm
Folded width	25 cm
Folded length	68 cm
Total weight	12 kg
Max. Load	130 kg
PU Pal.	8 Pal.