

Walking frame Grand XL

Item No. 10220HD

Image not found or t



- Rigid walking frame for maximum requirements
- The unique U-shape of the design enables greater mobility
- Wide, robust frame made of 2.5 cm anodized aluminum tube
- Offers maximum stability and light weight at the same time
- Can be used indoors and outdoors
- Space-saving folding by releasing the two locks

Technical Data

Order number	10220HD
Colour	white
Overall width	62 cm
Overall length	52.5 cm
Overall height	82 - 100 cm
Width between handles	52,5 cm

Folded width	62 cm
Folded height	82 cm
Folded depth	19 cm
Total weight	3.5 kg
Max. Load	228 kg
PU pc.	2 Stk.