

Oxygen concentrator 1025KS 10-liter

Item No. 1025KS



- Stationary 10-liter oxygen concentrator
- Oxygen concentration up to 96 % in the 2 - 10 l/min range
- Integrated OSD® sensor for continuous oxygen concentration measurement
- Visual and acoustic alarm for insufficient oxygen concentration, power failure, blocked oxygen outlet or device malfunction
- Front label with generally understandable image symbols and LED lamps
- Additional oxygen outlet on the rear side of the device; enables the use of all filling stations commonly available on the market in combination with an optional storage compartment
- Oxygen outlet equipped with fire protection adapter
- Multiple test runs completed successfully at very high temperatures and high humidity
- Manufactured in the USA

Technical Data

Order number	1025KS
Overall width	34.2 cm
Overall height	62.2 cm
Overall depth	30.4 cm
Flow rate	2 - 10 l/min

Oxygen concentration	87% - 96%
Power Supply	50 Hz
Total weight	19 kg
PU Pal.	18 Pal.

Oxygen concentrator 1025KS 10-liter

Item No. 1025KS



Accessories



Firesafe™ MKII cannula valve

Item No.827-2001x20