

Folding safety grab bar with fall protection

Item No. 10495-30

drive

DeVilbiss
HEALTHCARE



- Made of mirror polished stainless steel
- For fastening on both sides of the toilet
- Can be folded up
- Recessed grip grooves ensure a secure hold
- Including fall protection and tile protection stop
- Surface guarantees maximum possible freedom from germs and bacteria
- Mounting plate (LxW) 27.5 x 10 cm
- Delivery including mounting materials

Technical Data

Order number	10495-30
Overall length	60 cm
Mounting plate width	10 cm
Mounting plate length	27 cm

Tube diameter	3,2 cm
Total weight	2.9 kg
Max. Load	120 kg