

Portable oxygen concentrator iGo2

Item No. 125K



Through its groundbreaking SmartDose technology, the iGo2 has created a completely new class within portable oxygen concentrators. The iGo2 automatically adjusts the oxygen delivery bolus rate to the patient's changing breathing frequency.

- Patented SmartDose™ technology
- Extremely sensitive respiration detection (trigger)
- Max Oxygen production: 1100ml
- Approved for use on aircraft
- Including practical carrying bag
- Small, light, and handy

Technical Data

Order number	125K
HMV-Nr.	14.24.06.1011
Colour	gray
Overall width	8.9 cm
Overall height	21.3 cm
Overall depth	21.8 cm

Oxygen concentration	90+4/-3
Power Supply	50/60 Hz
Sound level	37 dBa
Total weight	2.25 kg
Weight	2.2 kg
PU Pal.	48 Pal.

Portable oxygen concentrator iGo2

Item No. 125K



Accessories



Battery charger for iGo2

Item No.125-CH-614



Rechargeable battery for iGo2

Item No.125D-613