

# Portable oxygen concentrator iGo2

Item No. 125K-alt



- Patented SmartDose™ technology
- Extremely sensitive respiration detection (trigger)
- Max Oxygen production: 1014ml
- Approved for use on aircraft
- Including practical carrying bag
- Small, light, and handy

The pioneering iGo2 has established an entirely new class of portable oxygen concentrators. It is the first of its kind to be equipped with the patented, self-adjusting SmartDose™ technology.

Oxygen output is automatically adapted to your activities. You can finally enjoy your life to the full again without having to constantly pay attention to the oxygen setting you need.

The patented SmartDose™ technology automatically adjusts the oxygen dose to your breathing pattern. More oxygen is output at times of increased activity.

The dose is reduced accordingly during your rest phases.

SmartDose™ therefore fosters an active and individual lifestyle.

For further information, please visit our website: [igo2poc.com](http://igo2poc.com)

## Technical Data

Order number	125K-alt
HMV-Nr.	14.24.06.1011
Colour	gray
Overall width	8.9 cm
Overall height	21.3 cm
Overall depth	21.8 cm

Oxygen concentration	90+4/-3
Power Supply	50/60 Hz
Sound level	37 dBa
Total weight	2.25 kg
Weight	2.2 kg
PU Pal.	48 Pal.

# Portable oxygen concentrator iGo2

Item No. 125K-alt

**drive**

**DeVilbiss**  
HEALTHCARE

## Accessories



**Battery charger for iGo2**

Item No.125-CH-614



**Rechargeable battery for iGo2**

Item No.125D-613