

Sling Split Leg in Chair Hammock

Item No. 4003XBE

drive

DeVilbiss
HEALTHCARE



- Ideal for decubitus and pain patients or for long periods of sitting
- Optimal patient comfort thanks to seamless bottom area and thinner quilting of the leg strap
- Reduces pressure and shear force thanks to the three-dimensional, multi-elastic 'spacer' knitted fabric
- Additional pressure relief through fold-away hip straps and thigh supports
- Max. load: 300 kg
- Machine-washable up to 80 °C

Technical Data

Order number	4003XBE
Size	XL
Total weight	1 kg

Weight	1 kg
Max. Load	300 kg