

Sling In Chair Hammock Comfort

Item No. 4101XBE



- Ideal for decubitus and pain patients or for long periods of sitting
- Optimal patient comfort thanks to seamless bottom area and thinner quilting of the leg strap
- Reduces pressure and shear force thanks to the three-dimensional, multi-elastic 'spacer' knitted fabric
- Max. load: 300 kg
- Machine-washable up to 80 °C

Technical Data

Order number	4101XBE
Size	XS
Total weight	1 kg

Weight	1 kg
Max. Load	300 kg