

Bath lift Bellavita 2G

Item No. 465200312

drive

DeVilbiss
HEALTHCARE



- New, unique 3D suction cups
- Service-friendly disassembly
- Maximum stability at highest position
- Wider side flaps available as accessories
- The lightest bath lift with a weight of just 9.3 kg
- Ergonomically shaped seat and back section
- Charge status is continuously indicated
- Front suction cups ventilate automatically on raising
- Tested for electromagnetic compatibility (EMC)
- No height adapters required (for bathtubs up to a depth of 48 cm)
- Five-year warranty for the bath lift, two-year warranty for the manual control and the battery (not including suction cups, covers, power supply unit, and parts subject to wear)

Article numbers in the European market

Please note the different article numbers valid outside Germany:

475200100 without cover

475200252 with cover Classic white

475200312 with cover Classic blue

475250252 with cover Comfort white

475250312 with cover Comfort blue

Technical Data

Bath lift Bellavita 2G

Item No. 465200312



Order number	465200312
Colour	blue
Overall width	70 cm
Overall length	83 cm
Seat width	38 cm
Seat depth	42 cm
Seat height	6.5 - 48 cm
Seat plate overall length	57 cm
Backrest width	32 - 36 cm
Backrest height	64 cm

Backrest angle	10 - 50 °
Baseplate width	31 cm
Baseplate length	56 cm
Total weight	9.7 kg
Weight	9.3 kg
Weight of battery-operated hand control	0,4 g
Seat unit weight	5,5 kg
Weight of backrest unit incl engine without battery	3,8 kg
Max. Load	140 kg
PU Pal.	30 Pal.

Bath lift Bellavita 2G

Item No. 465200312

drive

DeVilbiss
HEALTHCARE

Accessories



Swivel aid Vitaturn K

Item No.412100100



Swivel and transfer aid Vitaturn L

Item No.412100200



Active foam cleaner

Item No.410900100



Suction cup set 4 x gray

Item No.460900220



Suction cup set 4 x transparent

Item No.460900225



Side flaps, Bellavita widening

Item No.460900227



Side flaps, Bellavita widening

Item No.460900228



Side flaps, Bellavita widening

Item No.460900229



Cover set Bellavita

Item No.460800252



Cover set Bellavita

Item No.460800312



Cover set Bellavita

Item No.460900252



Cover set Bellavita

Item No.460900312



Side flap deflectors

Item No.400570



Belt systems, 4-point chest strap

Item No.460900150

Bath lift Bellavita 2G

Item No. 465200312



Belt systems, 2-point chest strap

Item No.460900100