

Bath lift Bellavita NOVA

Item No. 465400312

drive

DeVilbiss
HEALTHCARE



- Patented, detachable backrest, easy to remove
- Easy care seat with practical hygiene cutout
- Twice as safe: more stability thanks to two racks
- Stable: patented 3D suction cups
- Ergonomic backrest with inclination up to 50° degrees

First Intallation

Article numbers in the European market

Please note the different article numbers valid outside Germany:

475400252 with Classic white cover

475400312 with cover Classic blue

475500252 with cover Comfort white

475500312 with cover Comfort blue

Technical Data

Order number	465400312
Colour	white
Overall width	70 cm
Overall length	102 cm
Seat width	38 cm
Seat depth	48 cm
Seat height	6.5 - 48 cm
Seat plate overall length	60 cm
Backrest width	32 - 35 cm

Backrest height	64 cm
Backrest angle	10 - 50 °
Baseplate width	31 cm
Baseplate length	61 cm
Weight	9.3 kg
Weight of battery-operated hand control	0,4 g
Max. Load	140 kg
PU Pal.	30 Pal.

Bath lift Bellavita NOVA

Item No. 465400312

drive

DeVilbiss
HEALTHCARE

Accessories



Swivel aid Vitaturn K

Item No.412100100



Swivel and transfer aid Vitaturn L

Item No.412100200



Active foam cleaner

Item No.410900100



Suction cup set 4 x gray

Item No.460900220



Suction cup set 4 x transparent

Item No.460900225



Side flaps, Bellavita widening

Item No.460900227



Side flaps, Bellavita widening

Item No.460900228



Side flaps, Bellavita widening

Item No.460900229



Cover set Bellavita

Item No.460800252



Cover set Bellavita

Item No.460800312



Cover set Bellavita

Item No.460900252



Cover set Bellavita

Item No.460900312



Side flap deflectors

Item No.400570



Belt systems, 4-point chest strap

Item No.460900150

Bath lift Bellavita NOVA

Item No. 465400312



Belt systems, 2-point chest strap

Item No.460900100