

Raised toilet seat TSE 150

Item No. 5200500

drive

DeVilbiss
HEALTHCARE



- Two hygienic cutaways and integrated seat
- Recessed grip on the lid makes it easy to open
- Folding, non-slip armrests
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces simplify cleaning
- Delivery includes assembly gloves
- Delivery in practical take-away packaging

Technical Data

Order number	5200500
Colour	white
Height adjustment	6.5 / 10 / 13.5 cm
Overall width	46 cm
Overall length	48.5 cm
Overall height	22.5 cm

Opening width	23 cm
Opening length	25 cm
Total weight	3.7 kg
Max. Load	150 kg
PU Pal.	32 Pal.
PU pc.	4 Stk.