

# Raised toilet seat TSE 150 Plus

Item No. 520150000



- NEW: innovative top assembly
- Two hygienic cutaways and integrated seat
- Recessed grip on the lid makes it easy to open
- Folding, non-slip armrests
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces simplify cleaning
- Delivery includes assembly gloves
- Delivery in practical take-away packaging

## Technical Data

|                   |                    |
|-------------------|--------------------|
| Order number      | 520150000          |
| Colour            | white              |
| Height adjustment | 6.5 / 10 / 13.5 cm |
| Overall width     | 56 cm              |
| Overall length    | 48.5 cm            |
| Overall height    | 41 cm              |
| Opening width     | 23 cm              |

|                       |         |
|-----------------------|---------|
| Opening length        | 25 cm   |
| Width between armrest | 48 cm   |
| Total weight          | 5.7 kg  |
| Max. Load             | 150 kg  |
| PU Pal.               | 32 Pal. |
| PU pc.                | 4 Stk.  |