

Raised toilet seat TSE 120 Plus

Item No. 520500500



- NEW: with innovative top mounting
- Two hygienic cutaways and integrated seat
- Recessed grip on the lid makes it easy to open
- Folding, non-slip armrests
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces simplify cleaning
- Delivery includes assembly gloves
- Delivery in practical take-away packaging

Technical Data

Order number	520500500
Colour	white
Height adjustment	8.5 / 12 / 15.5 cm
Overall width	55 cm
Overall length	49 cm
Overall height	38 cm

Opening width	20,5 cm
Opening length	27 cm
Width between armrest	46 cm
Total weight	4.1 kg
Max. Load	120 kg
PU Pal.	30 Pal.