

Raised toilet seat Ticco 2G/5

Item No. 523600100

drive

DeVilbiss
HEALTHCARE



- Raises the toilet seat by 5 cm
- Ergonomic shape
- Maximum opening with two hygienic cutaways
- Recessed grip on the lid makes it easy to open
- Simple, tool-free assembly
- Fixation clamps enable perfect hold
- Smooth surfaces for optimal cleaning and disinfection
- Lightweight with maximum load capacity
- Delivery in practical take-away packaging

Technical Data

Order number	523600100
Colour	white
Height front	4,1 cm
Height rear	5,6 cm
Overall width	39 cm
Overall length	39 cm
Opening width	20,3 cm

Opening length	27,3 cm
Total weight	820 g
Weight	820 g
Max. Load	225 kg
PU Pal.	120 Pal.
PU pc.	12 Stk.