

Raised toilet seat Ticco 2G/15

Item No. 523600300

drive

DeVilbiss
HEALTHCARE



- Raises the toilet seat by 15 cm
- Ergonomic shape
- Maximum opening with two hygienic cutaways
- Recessed grip on the lid makes it easy to open
- Simple, tool-free assembly
- Fixation clamps enable perfect hold
- Smooth surfaces for optimal cleaning and disinfection
- Lightweight with maximum load capacity
- Delivery in practical take-away packaging

Technical Data

Order number	523600300
Height front	13,6 cm
Height rear	15,1 cm
Overall width	39 cm
Overall length	39 cm
Opening width	20,3 cm

Opening length	27,3 cm
Total weight	1.44 kg
Max. Load	225 kg
PU Pal.	60 Pal.
PU pc.	6 Stk.