

Oxygen concentrator Compact 525 5-liter

Item No. 525KS



- Stationary 5-liter oxygen concentrator
- Increases patient comfort while simultaneously reducing operating costs
- Oxygen concentration 93 % \pm 3 % at 0.5 - 5 l/min
- Integrated OSD® sensor for continuous oxygen concentration measurement
- Visual and acoustic alarm for insufficient oxygen concentration, power failure, blocked oxygen outlet or device malfunction
- Intelligent Turn-Down® technology reduces the concentrator's power consumption
- The 525KS is very easy to move thanks to the attached casters and its integrated handle
- Manufactured in the USA

Technical Data

Order number	525KS
Colour	black
Overall width	34.2 cm
Overall height	62.2 cm
Overall depth	30.4 cm

Flow rate	0.5 - 5 l/min
Oxygen concentration	93 % \pm 3 %
Power Supply	50/60 Hz
Total weight	16.3 kg
PU Pal.	18 Pal.

Oxygen concentrator Compact 525 5-liter

Item No. 525KS

Accessories



O2 Starter Set

Item No. O2-Starterset



Firesafe™ MKII cannula valve

Item No. 827-2001x20