

Wheeled commode TRS 130

Item No. 530100000

drive

DeVilbiss
HEALTHCARE



- Comfortable, safe, and particularly stable
- Parking brake on the rear casters
- Fold-down armrests simplify transfer
- Stable, folding, and removable footrests
- Receptacle with odor trap as standard
- Can be positioned over the toilet bowl without receptacle
- Continuously adjustable footrests
- Simple, tool-free assembly
- Calf strap. Easy to clean

Technical Data

Order number	530100000
Overall width	56 cm
Overall length	95 cm
Overall height	97 cm
Seat width	44 cm
Seat depth	42.5 cm
Seat height	- 52 cm

Backrest height	38 cm
Clearance height	43 cm
Footrest adjustment	42 - 50 cm
Total weight	14.4 kg
Weight	14.4 kg
Max. Load	130 kg
PU Pal.	16 Pal.

Wheeled commode TRS 130

Item No. 530100000



Accessories



Bucket with lid

Item No.530100090



Seat for TRS 130

Item No.530100001



Seat cover for TRS 130

Item No.530100020