

Wheeled commode TRS 200 XXL

Item No. 530200000

drive

DeVilbiss
HEALTHCARE



- Wide and stable wheeled commode for maximum requirements
- Extra-large and soft seat surface and backrest for increased comfort
- Seat and cover can be folded up
- Parking brake on the rear casters
- Fold-down armrests simplify transfer
- Stable, folding, and removable footrests
- Continuously adjustable footrests
- Extremely simple assembly and cleaning
- With receptacle as standard

Technical Data

Order number	530200000
Overall width	72 cm
Overall length	96 cm
Overall height	98 cm
Seat width	56 cm
Seat depth	42 cm
Seat height	- 52 cm

Clearance height	44 cm
Footrest adjustment	42 - 57 cm
Total weight	15.3 kg
Weight	15.3 kg
Max. Load	200 kg
PU Pal.	6 Pal.

Wheeled commode TRS 200 XXL

Item No. 530200000

Accessories



Seat for TRS 200 XXL

Item No.530200021



Seat cover for TRS 200 XXL

Item No.530200026