

Oxygen filling station iFill

Item No. 5351



- The iFill oxygen filling station ensures an unlimited supply of refillable oxygen
- The user-friendly device is simple to use and operates reliably
- The filling station can be positioned anywhere in the home
- Various bottle sizes available for all requirements
- Bottles can be filled as often as required by the user
- No more waiting for your bottles to be delivered
- iFill kit consisting of iFill, power cord, 2 oxygen bottles (type depending on configuration), a carrying bag, and instructions for use

Technical Data

Order number	5351
Colour	black
Overall width	57 cm
Overall height	72.5 cm

Durchschnittlicher Fülldruck	138 ± 13.8 bar 2000 ± 200 psi bar
Oxygen concentration	93 % ± 3 %
Power Supply	50/60 Hz
Total weight	30 kg