

Shower stool Delphi

Item No. 540100000

drive

DeVilbiss
HEALTHCARE



- Safe, stable, and space-saving
- Seat surface with drainage slits
- Low weight
- Three-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean
- Delivery in practical take-away packaging

Technical Data

Order number	540100000
Colour	white
Overall width	48 cm
Overall length	48 cm
Overall height	55 cm
Seat plate diameter	38 cm

Seat height	45 - 55 cm
Seat height adjustment	3
Total weight	1.6 kg
Max. Load	180 kg
PU Pal.	80 Pal.
PU pc.	10 Stk.