

Shower stool Delphi 2 Plus

Item No. 540300000



- With 2 hygiene openings for comfortable body care
- Safe, stable, and space-saving
- Seat surface with drainage slits
- Low weight
- Three-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean
- Delivery in practical take-away packaging

Technical Data

Order number	540300000
Colour	white
Overall width	48 cm
Overall length	48 cm
Overall height	55 cm
Seat plate diameter	38 cm
Seat height	45 - 55 cm

Seat height adjustment	3
Total weight	1.6 kg
Weight	1.6 kg
Max. Load	180 kg
PU Pal.	80 Pal.
PU pc.	10 Stk.