

Shower chair DSF 130

Item No. 550200100

drive

DeVilbiss
HEALTHCARE



- Space-saving thanks to folding mechanism
- With armrest and backrest
- Large seat surface with drainage holes
- Hygienic cutaway simplifies personal hygiene
- Comfortable, safe, and stable
- Five-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Powder-coated aluminum frame
- Easy to clean

Technical Data

Order number	550200100
Colour	white
Overall width	52 - 53 cm
Overall height	90 - 100 cm
Overall depth	56 - 58 cm
Seat width	41 cm
Seat depth	46 - 41 cm

Seat height	42 - 48 cm
Seat height adjustment	5
Backrest height	48 cm
Total weight	3.8 kg
Max. Load	130 kg
PU Pal.	20 Pal.
PU pc.	2 Stk.