

Shower stool Dino

Item No. 550500100

drive

DeVilbiss
HEALTHCARE



- Safe, stable, and space-saving
- Simple, tool-free assembly
- Very light
- Structured surface ensures secure hold
- Eight-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean

Technical Data

Order number	550500100
Colour	white
Overall width	45 cm
Overall length	48 cm
Overall height	52 cm
Seat plate diameter	32 cm

Seat height	34 - 52 cm
Seat height adjustment	8
Total weight	1.6 kg
Max. Load	130 kg
PU Pal.	120 Pal.
PU pc.	6 Stk.