

Swiveling bath seat BEN 2G

Item No. 550600300

drive

DeVilbiss
HEALTHCARE



- The seat can be swiveled to the desired position in 90° steps
- Automatic locking at 90° ensures safety
- Fully assembled
- Simply hook in at the edge of the bathtub
- Simple adjustment to the bathtub width (max. inner bathtub width 66 cm)
- Drainage holes allow the shower water to flow away
- Mounting arms with non-slip rubber caps
- Powder-coated steel frame with stainless steel pivot pins
- Easy to clean

Technical Data

Order number	550600300
Colour	white
Overall width	71 cm
Overall depth	52,8 cm
Seat width	44 cm
Seat depth	44 cm
Seat adjustment	4 x 90°

Adjustment (width of bathtub)	48.5 - 64.5 cm
Backrest height	37 cm
Total weight	3.47 kg
Max. Load	130 kg
PU Pal.	16 Pal.
PU pc.	2 Stk.