

Shower chair DH 135

Item No. 550800100

drive

DeVilbiss
HEALTHCARE



- Available with and without backrest
- Large seat surface with maximum possible hygienic cutaway
- Easy to transport and stow
- Stable seat in the bath and under the shower
- Comfortable, safe, and stable
- Simple design
- Can be dismantled
- Five-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean

Technical Data

Order number	550800100
Colour	white
Overall width	42.5 cm
Overall height	43 - 53 cm
Overall depth	47 cm
Seat width	41 cm
Seat depth	37 cm

Seat height	43 - 53 cm
Seat height adjustment	5
Total weight	2 kg
Max. Load	135 kg
PU Pal.	45 Pal.
PU pc.	5 Stk.