

Shower Chair AiO

Item No. 550900400

drive

DeVilbiss
HEALTHCARE



The modular AiO shower chair can be individually adapted to your needs and can be extended at any time as your life changes. It is characterized by its modular plug-in system and is therefore the perfect solution for a wide range of requirements. If your circumstances change, the AiO can be easily extended with additional components such as a high backrest or sturdy armrests. In no time at all, the AiO transforms from a simple shower stool into a comfortable shower chair that meets all your needs.

- Individually adaptable to your needs
- Beading effect simplifies care and cleaning
- Super lightweight yet robust
- Made of high-quality, skin-friendly plastic, easy to clean and corrosion-free
- Adjustable to your size thanks to 8-way adjustable seat height
- Non-slip rubber feet ensure stability
- Large hygiene cut-out for easy personal hygiene

Technical Data

Order number	550900400
Colour	white
Overall width	51,5 cm
Overall height	77,5 - 94,5 cm
Overall depth	46 - 48 cm
Seat width	42,7 cm
Seat depth	34 cm

Seat height	38 - 55 cm
Seat height adjustment	8
Backrest height	41 cm
Total weight	4.44 kg
Max. Load	130 kg
PU Pal.	30 Pal.
PU pc.	1 Stk.