

Child bath support Otter

Item No. 555100100

drive

DeVilbiss
HEALTHCARE



- Plastic frame with removable, water-permeable vinyl cover
- Seat surface and backrest can be adjusted to four positions (inclination from 0° to 90°)
- Safety straps for body and legs
- Head/neck support height- and width-adjustable; also suitable as side protection
- Can be stowed compactly when folded
- With stand or rolling frame, extensible for many everyday situations
- Head support, chest strap, and 2 leg straps are included in delivery

Technical Data

Order number	555100100
Colour	white-green
Body size	bis 90 cm
Seat width	41 cm
Seat depth	30 cm
Seat height	14 - 26 cm
	0 - 90 °
Backrest width	41 cm

Backrest height	40 cm
Backrest angle	0 - 90 °
Folded height	15 cm
Folded width	42,5 cm
Folded length	90 cm
Weight	5.1 kg
Max. Load	30 kg
PU Pal.	20 Pal.