

Shower support base frame, movable

Item No. 555200200

drive

DeVilbiss
HEALTHCARE



- Quickly transforms the Otter bath support into a shower support
- Base frame with casters for increased flexibility
- Pleasant height for the supervising person
- Suitable for all three sizes
- Robust aluminum frame
- Parking brake for the casters
- Click system enables easy set-up and dismantling

Technical Data

Order number	555200200
Colour	white-green
Overall width	69 cm
Overall length	95 cm

Overall height	52,5 cm
Weight	7 kg
Max. Load	85 kg