

Universal sling comfort with head support

Item No. 610000023



- The universal sling
- Quick to put on
- Ideal for patients who require head support
- Good support in the back and pelvis area
- Made of polyester
- Max. load: 200 kg
- Machine-washable up to 80 °C

Slings are available in various materials for different patients and environments:

Mesh - soft and pliable material, used primarily for bathing as the material allows water to drain away and dries quickly after use.

Polyester - a popular and standard material for slings. Polyester is robust yet smooth; it has low friction which makes the sling easy to put on and take off.

Technical Data

| | |
|--------------|-----------|
| Order number | 610000023 |
|--------------|-----------|

| | |
|-----------|--------|
| Max. Load | 200 kg |
|-----------|--------|