

Full body sling comfort with head support

Item No. 610000042



- The Hammock full body sling offers optimal support in the sacral and thigh area
- Also suitable for leg amputees
- Quick to put on
- Ideal for patients who require head support
- Additional quilting at the bottom for high comfort
- Made of polyester
- Max. load: 200 kg

Slings are available in various materials for different patients and environments:

Mesh - soft and pliable material, used primarily for bathing as the material allows water to drain away and dries quickly after use.

Polyester - a popular and standard material for slings. Polyester is robust yet smooth; it has low friction which makes the sling easy to put on and take off.

Technical Data

Order number	610000042
--------------	-----------

Max. Load	200 kg
-----------	--------