

Suction cup grab bar Solido T

Item No. 660100500

drive

DeVilbiss
HEALTHCARE



- The mobile, telescopic grab bar Solido T fits anywhere
- Makes it easier to hold on and stand up in the bathtub or shower
- Versatile and quickly ready to use on all smooth surfaces
- No screwing or drilling
- Vacuum rocker lever generates extra-strong suction with a tensile force of 70 kg
- High-quality finish and materials
- Articulated adapter available as an accessory
- Delivery in practical take-away packaging

Technical Data

Order number	660100500
Overall length	56 - 69 cm
	33 - 46 cm
Wall clearance	3,5 cm
Total weight	900 g

Weight	900 g
Max. Load	70 kg
PU Pal.	180 Pal.
PU pc.	10 Stk.

Suction cup grab bar Solido T

Item No. 660100500



Accessories



Articulated adapter for Solido T

Item No. 660100700